[www.zimed.com.tr](http://www.zimed.com.tr)

**zimed**<sup>®</sup>

**ALT-ÜST EKSTREMİTE KEMİK  
PLAKLARI KULLANIM KILAVUZU**  
*LOWER-UPPER EXTREMITY BONE  
PLATES INSTRUCTIONS FOR USE*

**TR**

**ENG**

**ES**

**IT**

**FR**

**RU**

*innovation to move*

# ZİMED MEDİKAL

## ALT-ÜST EKSTREMİTE KEMİK PLAKLARI KULLANIM KILAVUZU

### Uzman Hekimin İlgisine (TR)

Doküman No: IFU.01.01/00  
14.06.2022

#### (TR) ÜRÜN TANIMI:

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#### (EN) PRODUCT DEFINITION:

It includes plates, screws, and hand tools that will enable the application of these products to ensure the anatomical alignment of the broken bones, to keep them together, and to provide splinting with the ability to ensure the early mobility of the patient. The products are non-sterile and the implants in the set are for single use. See product catalog for reference numbers and names of products.

Lower-Upper Extremity Bone Plates are made of titanium alloy (ISO 5832-2, ISO 5832-3, ASTM F 67, ASTM F 136).

#### SMALL FRAGMENT PLATE SYSTEM;

**Locked LC DCP Narrow Plate 3.5 mm:** It is designed for superior flexing endurance, strong fracture determination. Thanks to its low-profile design, it also eliminates the problem of tissue irritation. It is indicated for fracture fixation and post-fixation osteotomy, malunion, nonunion cases, radius and ulna diaphysis (shaft) fractures. It is manufactured from titanium alloy produced according to Locked LC DCP Narrow Plate ASTM F136. It is used with Ø3.5 mm locked, Ø3.5 Cortical Screw. 4-14 holes and 72, 85.5, 99, 114, 125, 144, 159, 174.5, 190, 205, 220 mm length options are available.

**Locked LC DCP Large Plate 3.5 mm:** The Locked LC DCP Plate (Limited Contact Dynamic Compression Plate) System is designed to be used in the upper extremity region, especially for the fixation of small bones such as the distal humerus and radius. Ø3.5 mm cortical screw is used together with Ø 3.5 mm locked screw. Locked LC DCP Broad is produced from titanium alloy manufactured according to ASTM F136. 4-14 holes, 72, 85.5, 99, 112.5, 129, 144, 159, 174, 190, 205, 220 mm length options are available.

**Locked 1/3 Tubular Plate:** 1/3 Tubular Plates are a compression-capable titanium plate system with low-contact surface. These plates, which have been given anatomical features according to the region to be used, prevent soft tissue irritation and provide maximum fracture determination. It is indicated for ulna and radius diaphysis (shaft) fractures, fibula, distal and diaphysis fractures, olecranon and open fractures, displaced fractures, transversal fractures. 4-14 holes and 51, 63, 75, 87, 99, 111, 123, 135, 147, 159, 171 mm length options are available. Locked 1/3 Tubular Plate is produced according to ISO 5832-2 TiGr3 (ASTM F 67) standard. It is used with Ø3.5 mm cortical screw, Ø 3.5 mm locked screw.

**Locked Distal Radius Volar Plate:** Usage of Distal Radius is indicated for intra- or extra-articular fractures accompanied by possible metaphyseal defect, open joint reconstructions, fractures of the carpal and metacarpal region with distal radius, and corrective osteotomies. Distal Radius Locked Volar Plates are used in the non-unions, osteotomies and surgical treatment of fractures of the distal radius. Volar plates are used with Ø3.5 mm cortical screws and Ø2.4 mm. 12, 13, 14, 15, 16, 17 holes options are available. It is designed differently as right and left. Locked Distal Radius Volar Plates are produced from the material of ISO 5832-3 Ti6Al4V ELI (ASTM F 136) material.

**Locked Distal Radius Volar Extension Plate:** Zimed Medikal Distal Radius Plate extension piece allows the extension the length of the plates by the surgeons to treat the fracture without using another plate. The extension piece is designed in three different lengths as 59,2 mm, 89 mm, 119 mm, with 6, 9, 12 holes. There is one connecting screw to ensure the connection between the extension and the plate. Volar extension plate is available as right and left options. It is used with Ø3.5 mm cortical screw and Ø2.4 mm.

**Locked Distal Ulna Plate:** It is indicated for the surgical treatment of distal ulna fractures resulting in loss of stabilization of the distal radioulnar joint, distal ulna head fractures in which the articular surface is displaced, rotated or angled, and comminuted ulna neck fractures that threaten distal radioulnar joint stabilization. It is used with a locked distal ulna plate, Ø3.5 mm locked screw, Ø3.5 cortical screw, Ø2.4 mm locked screw. 7, 8, 9 holes, 46,6, 60, 68 mm length options are available. It has right and left options. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Distal Lateral Humeral Plate:** It is indicated for the surgical treatment of non-unions of the anatomical region associated with intra-articular and supracondylar fractures of the distal humerus. Locked Distal Lateral Humeral Plate is used with Ø3.5 mm locked screw, Ø3.5 mm Cortical screw, Ø4.0 mm locked spongiosis screw, Ø4.0 locked cannulated

spongiosis screw. 8, 10, 12, 14 holes, 86, 106, 126, 146 mm length options are available. It is designed differently as right and left. It is produced from ISO 5832-3 Ti6Al4V ELI (ASTM F 136) material.

**Locked Distal Medial Humeral Plate:** It is indicated for the fracture fixation and post-fixation osteotomy of the distal humerus, non-union surgical treatment. 8, 10, 12, 14, 16 holes, 87, 107, 127, 147, 167 mm length options are available. Right-left feature is available. Locked Distal Medial Plate is produced from ISO 5832-3 Ti6Al4V ELI (ASTM F 136) material. It is used with Ø3.5 mm locked screw, Ø3.5 mm Cortical screw, Ø4.0 mm locked spongiosis screw.

**Locked Distal Radius Plate Dorsal:** Distal Radius is indicated for extra-articular fractures with possible metaphyseal defect, open joint reconstructions, fractures of the carpal and metacarpal region seen with distal radius, and corrective osteotomies. Distal Radius Locked Plates are used in the distal radius non-unions and surgical treatment of fractures. It is designed as right and left. It is used with Ø3.5 mm cortical screw, Ø 3.5 mm locked screw and 2.4 locked screw. 10-15 holes, 51, 59, 67, 80, 88, 96 mm length options are available. Locked Distal Radius Dorsal Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Distal Posterior Humeral Plate:** It is indicated for the fixation of distal posterior humeral joint fractures and fixation of osteotomies. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material. 8-12 holes, 97, 110, 123, 136, 150 mm length options are available. It is used with Ø 3.5 mm locked screw and Ø 3.5 mm Cortical screw. It has right and left options.

**Locked Proximal Lateral Humeral Plate (Large):** It is indicated for the surgical treatment of non-unions of the anatomical region associated with intra-articular and supracondylar fractures of the proximal humerus. It is produced from ISO 5832-3 Ti6Al4V ELI (ASTM F 136) material. 12, 14, 16, 18, 20 holes, 104, 133.5, 163, 193.5, 222.5 mm length options are available. It is used with Ø3.5 locked screw, Ø3.5 cortical screw, Ø3.5 locked cannulated screw, Ø4.0 locked spongiosis screw, Ø4.0 locked cannulated spongiosis screw.

**Locked Proximal Radial Head Plate:** The Radial Head System is indicated for proximal radius fractures, fusions and osteotomies. Used with Ø3.5 mm locked screw, Ø2.4 locked screw, Ø3.5 cortical screw. 6, 7, 8 holes 34, 42, 50 mm length options are available. Locked Proximal Radial Head Plate is manufactured from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Distal Humeral Y Plate:** It is indicated for intra-articular fractures of the distal humerus, supracondylar fractures, and non-union conditions of the distal humerus. It is used with Ø3.5 locked screw, Ø3.5 cortical screw. Locked Distal Humeral Y Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material. 9-12 holes, 73, 85, 97, 109, 121, 133, 145, 157 mm length options are available. It is designed as right and left.

**Locked Wrist Arthrodesis Plate:** The Locked Wrist plate system is indicated for wrist arthrodesis and fractures of other small bones of the carpus. It is used with Ø3.5 mm locked screw and Ø3.5 cortical screw. 8 holes size option are available. Locked Wrist Arthrodesis Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Distal Clavicle Plate:** The Distal Clavicle Plate is indicated for fractures close to the distal portion of the Clavicle bone and for fractures in the distal portion and osteotomies. It is used with Ø 2.4 mm locked screw, Ø3.5 mm locked screw and Ø3.5 cortical screw. 14 holes and 17 holes right and left size options are available. Distal Clavicle Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Clavicle Plate:** The Locked Clavicle Plate is indicated for fractures and osteotomies of the shaft part of the Clavicle bone. It is used with Ø3.5 mm locked screw and Ø3.5 cortical screw. 6, 8, 10, 12 holes right and left size options are available. Clavicle Plate is produced from Ti6Al4V ELI (ASTM F 136) material.

**Locked Olecranon Plate:** The Locked Olecranon Plate is indicated for complex extra-articular and intra-articular olecranon fractures, pseudoarthrosis of the proximal ulna. It is used with Ø3.5 mm locked screw and Ø3.5 cortical screw. 8, 10, 12, 14, 16 holes right and left size options are available. Olecranon Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Straight Reconstruction Plate 3.5 mm:** Locked Straight Acetabulum Reconstruction Plate is designed to be used to fix bones in fractures of the acetabulum region. It is used with Ø3.5 mm cortical screw, Ø 3.5 mm locked screw, Ø 3.5 mm locked cannula and 4.0 mm spongiosis screw. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material. There are different length options as 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20 holes.

**Locked Clavicle Hook Plate:** The Locked Clavicle Hook Plate is indicated for dislocation of the acromioclavicular joint and fixation of lateral Clavicle fractures. It is used with Ø3.5 mm locked screw and Ø3.5 cortical screw. 5, 7, 9 holes right and left size options are available. Clavicle Hook Plate, Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Distal Fibula Plate:** The Locked Distal Fibula Plate is indicated for fractures and osteotomies of the distal fibula metaphyseal region of the Fibula bone. It is used with Ø3.5 mm locked screw and Ø3.5 cortical screw. 9, 11, 13 holes right and left size options are available. Distal fibula plates are manufactured from titanium alloy produced according to ASTM F 136.

**Locked Distal Lateral Fibula Plate:** The Locked Distal Lateral Fibula Plate is indicated for fractures of the distal lateral fibular metaphyseal region of the Fibula bone and osteotomies. It is used with Ø2.4 mm locked screw, Ø3.5 mm locked screw and Ø3.5 cortical screw. 9, 11, 13, 15, 17 holes right and left size options are available. Distal lateral fibula plates are manufactured from titanium alloy produced according to ASTM F 136.

**Eight Growth Plates:** The Eight Growth Plates are indicated for correcting angular growth deformities in children, Femur and/or tibia extension deformities, and other medium and small bone deformities. It is used with Ø3.5 cortical screw. 12 mm, 16 mm, 20 mm and 24 mm length options are available. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Calcaneus Plates:** It is indicated for fixation of complicated calcaneus fractures, calcaneus comminuted fractures and osteotomies. Locked lateral calcaneus plate is available with options 78 mm and 95 mm length and 18 and 20 holes; locked calcaneus Y plate with options 46 mm, 60 mm, 70 mm length and 8 holes; calcaneus F plate with options 46 mm, 60 mm, 70 mm length and 11, 14 holes; right and left options are available. Calcaneus plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material and used with Ø 4.0 mm locked spongiosis, Ø 4.0 mm spongiosis screw and Ø3.5 locked and cortical screw.

**Locked Proximal Lateral Femoral Plate 3.5 mm:** The Proximal Lateral Femoral Plate is indicated for trochanteric region fractures, trochanterodiaphyseal, proximal and shaft fractures seen together, malunion / nonunion conditions, and proximal femur osteotomies. It is used with Ø3.5 mm locked screw, Ø4.5 cortical screw and Ø6.5 spongiosis screw. 9, 11, 13, 15, 17 holes right and left size options are available. The plate is produced from Ti6Al4V ELI (ASTM F 136) material.

#### LARGE FRAGMENT PLAQUE SYSTEM;

**Locked LC DCP Broad Plate 4.5 mm:** It is designed for strong fracture fixation. Thanks to its low profile design, it also eliminates the problem of tissue irritation. It is indicated in the case of fracture fixation and post-fixation osteotomy, non-union, in diaphysis (shaft) fractures. It is manufactured from titanium alloy produced according to Locked LC DCP Broad Plate ASTM F 67 and used with Ø 5 mm Locked Screw and Ø 4,5 Cortical Screw. Size options are available between 4-16 holes.

**Locked LC DCP Narrow Plate 4.5 mm:** It is designed for strong fracture fixation. Thanks to its low profile design, it also eliminates the problem of tissue irritation. It is indicated for fracture fixation and post-fixation osteotomy, malunion, non-union cases, Tibia diaphysis (shaft) fractures. Locked LC DCP Narrow Plate is manufactured from titanium alloy produced

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Sti.

Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye

Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

# ZİMED MEDİKAL

## LOWER-UPPER EXTREMITY BONE PLATES INSTRUCTIONS FOR USE

### Submitted to the attention of Surgeons

### (EN)

Document No: IFU.01.01/00  
14.06.2022

according to ASTM F 67. It is used with Ø 5 mm Locked Screw and Ø 4.5 Cortical Screw. Size options are available between 4-16 holes.

**Locked Distal Lateral Tibial Plate:** Locked Distal Lateral Tibial Plate is indicated for extra-articular and simple intra-articular distal tibia fractures, distal tibial fractures extending to the diaphysis region. It is used with Ø6.5 mm locked spongiosis screw, Ø5 mm locked screw and Ø4.5 cortical screw. 8, 10, 12, 14, 16 holes right and left size options are available. Distal Lateral Tibia Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Proximal Medial Tibial Plate:** The Proximal Medial Tibial Plate is indicated for trochanteric region fractures, trochanterodiaphyseal, in fractures where proximal and shaft fractures seen together, malunion / nonunion conditions, and proximal tibial osteotomies. It is used with Ø5mm locked screw, Ø4.5 cortical screw and Ø6.5 spongiosis screw. 8, 10, 12, 14, 16 mm holes right and left size options are available. The plate is produced from titanium material manufactured according to ISO 5832-2.

**Locked Distal Medial Tibial End Point Plate (3.5-5.0 mm combined holes):** Locked Distal Medial Tibial End Point Plate is indicated for Complex extra-articular and intra-articular distal tibial fractures, and osteotomies. It is used with Ø5mm locked screw, Ø 3.5 locked screw, Ø 4.5 cortical screw, Ø5.0 mm locked screw and Ø4.0 locked or non-locked spongiosis screw. 13, 15, 17, 19, 21 holes right and left size options are available. Distal Medial Tibia End Point Plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Proximal Lateral Femoral Plate 5.0 mm:** Proximal Lateral Femoral Plate is indicated for trochanteric region fractures, trochanterodiaphyseal, proximal and shaft fractures seen together together, malunion / non-union conditions and proximal femur osteotomies. It is used with Ø5mm locked screw, Ø4.5 cortical screw and Ø6.5 spongiosis screw. 9, 11, 13, 15, 17 holes right and left size options are available. The plate is produced from Ti6Al4V ELI (ASTM F 136) material.

**Locked Proximal Lateral Femoral Head Plate:** It is used in the surgery of proximal radius head fractures, for the fixation of fractures close to the proximal radius, fractures with radius proximal comminuted fractures, in radius proximal corrective osteotomy. It is used with Ø4.5 mm cortical screw and Ø7.3 mm locked cannulated screw suitable for the anatomical structure of the plate radius head. It has right and left features with 4 holes options. It is produced from Ti6Al4V ELI (ASTM F 136) material.

**Locked Proximal Medial Tibial T Plate:** It is indicated for proximal tibial fractures. It is used with Ø6.5 mm spongiosis screw, Ø4.5 mm cortical screw and Ø5.0 mm locked screw. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Locked Proximal Lateral Tibial L Plate:** It is indicated for the surgical treatment of metaphyseal and split type fractures of the tibial plateau. It is used with Ø6.5 mm spongiosis screw, Ø4.5 mm cortical screw and Ø5.0 mm locked screw. The plate is produced from Ti6Al4V ELI (ASTM F 136) material.

**Locked Distal Lateral Femoral Plate:** The Distal Lateral Femoral Plate is indicated for Supracondylar fractures, intra-articular and extra-articular condylar fractures, defective union and non-union of the distal femur, and periprosthetic fractures. It is used with Ø5mm locked screw, Ø4.5 cortical screw and Ø6.5 Spongiosis screw. The plate has 13, 15, 17, 19, 21 holes right and left size options. The plate is produced from Ti6Al4V ELI (ASTM F 136) material.

**Adjustable Locked High Tibial Osteotomy Plate Type 1:** High Tibial Locked Osteotomy Plate is indicated for open wedge/closed wedge osteotomies of the medial proximal tibia, idiopathic or posttraumatic varus or valgus deformities of the proximal tibia. It is used with Ø5mm locked screw and Ø6.5 spongiosis screw. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Adjustable Locked High Tibial Osteotomy Plate Type 2:** It is used with Locked High Tibial Osteotomy Plate. There are 6, 8, 10, 12, 12 mm length options. It is produced from Ti6Al4V ELI (ASTM F 136) material.

#### Locked Distal Anterolateral Tibial Plate:

Locked distal Anterolateral tibial plate is indicated for fixation of distal tibial fractures, osteotomies and non-unions of patients with osteopenic bone tissue. It is used with Ø3.5 locked screw, Ø3.5 cortical screw, Ø4.0 spongiosis screw. It is produced from locked distal anterolateral tibial Ti6Al4V ELI (ASTM F 136) material. There are 10, 12, 14, 16, 18 holes and right and left options.

**Patella Plate:** It is indicated for the correction of patella fractures and osteotomies. S, M, L size options are available. Patella plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

#### MINIMAL INVASIVE LOCKED PLATES;

**Minimally Invasive Locked Proximal Lateral Humeral Plate (Small);** Minimally invasive proximal lateral humeral plate is used with guide system. It is applied by sending a screw to the plate over the guide with a very small incision compared to normal practice. It is indicated for stabilization of proximal humeral fractures. It is used with Ø3.5 mm locked and cortical screw, Ø4.0 mm locked spongiosis screw and Ø4.0 mm spongiosis screw. The plate is produced from Ti6Al4V material.

**Minimally Invasive Locked Distal Lateral Femoral Plate;** The minimally invasive Distal Lateral Femoral Plate is used with the guide system. It is applied by sending a screw to the plate over the guide with a very small incision compared to normal practice. It is indicated for supracondylar fractures, intra-articular and extra-articular condylar fractures, defective union and non-union of the distal femur, and periprosthetic fractures. It is used with Ø5mm locked screw. The plate has 12, 16, 20 holes right and left size options. The plate is produced from Ti6Al4V material.

**Minimally Invasive Locked Proximal Lateral Tibial Plate;** Minimally invasive Proximal Lateral Tibial Plate is used with guide system. It is applied by sending a screw to the plate over the guide with a very small incision compared to normal practice. It is indicated for stabilization of proximal tibia fractures (Proximal shaft fractures, Metaphyseal fractures, intra-articular fractures and periprosthetic fractures). It is used with Ø5mm self drilling locked screw. The plate is available in 10, 14, 18 holes right and left size options. The plate is produced from Ti6Al4V material.

#### ANGLE BLADE PLATES;

##### Angled Wedge Plate 90°:

The plate is indicated for stabilization and straightening of bone in the distal femur. It is used with Ø3.5mm locked screw and Ø3.5 cortical screw. 3, 4, 5, 6, 7, 8 holes length options and 28 mm, 33 mm, 38 mm wedge lengths are available. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

##### Locked Angle Wedge Plate 135°:

The plate is indicated for stabilization and straightening of bone in the proximal femur. It is used with Ø3.5mm locked screw and Ø3.5 cortical screw. 3, 4, 5, 6, 7, 8 holes length options and 28 mm, 33 mm, 38 mm wedge lengths are available. The plate is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

#### PELVIS PLATE SYSTEM

**Locked Straight Reconstruction Plate:** The Reconstruction Plate System is designed to be used in the fixation of bones in bone fractures of the Pelvis region. Used with Ø3.5 mm cortical screw. There are different size options between 3-20 holes. The plate is produced from titanium alloy manufactured according to ISO 5832-2.

**Curved Reconstruction Plate:** The Curved Reconstruction Plate System is designed to be used in the fixation of bones in pelvic region bone fractures. Used with Ø3.5 mm locked screw. There are different size options between 4-16 holes. The plate is produced from titanium alloy manufactured according to ISO 5832-2.

**Acetabular Hook Plate:** Pelvic System consists of plates for the treatment of sacroiliac joint dislocations and symphysis pubis disruptions during fractures, fusions and osteotomies of the acetabulum, sacrum, ilium and the entire pelvic ring. 2, 3, 4, 5, 6 holes options are

available in different sizes. It is used with Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Posterior Wall acetabular plate:** Pelvis System consists of plates for the treatment of sacroiliac joint dislocations and symphysis pubis disruptions during fractures, fusions and osteotomies of acetabulum, sacrum, ilium and the entire pelvic ring. 9, 10 holes right-left options are available. Used with Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Symphysis Pubis Plate,** Pelvis System consists of plates for the treatment of sacroiliac joint dislocations and disruptions of the symphysis pubis during fractures, fusions and osteotomies of acetabulum, sacrum, ilium and the entire pelvic ring. 4 holes S, M, L options are available. It is used with Ø 3.5 mm cortical Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Illipectineal Pelvis Plate,** Pelvis System consists of plates for the treatment of sacroiliac joint dislocations and symphysis pubis disruptions during fractures, fusions and osteotomies of acetabulum, sacrum, ilium and the entire pelvic ring. 12 holes right and left options are available. It is used with Ø 3.5 mm cortical Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Anterior Brim Pelvis Plate,** Pelvis System consists of plates for the treatment of sacroiliac joint dislocations and symphysis pubis disruptions during fractures, fusions and osteotomies of acetabulum, sacrum, ilium and the entire pelvic ring. 13 holes right and left options are available. It is used with Ø 3.5 mm cortical Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

**Superior Sacroiliac Plate,** Pelvis System consists of plates for the treatment of sacroiliac joint dislocations and symphysis pubis distortions during fractures, fusions and osteotomies of acetabulum, sacrum, ilium and the entire pelvic ring. 4 holes one size option is available. It is used with Ø3.5 mm cortical screw and Ø 3.5 mm locked screw. Plates are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

#### SCREWS USED:

**Ø5.0 mm Locked Screw:** Ø 5.0 mm locked screws are applied together with locked plates. It is used to fix the lower extremity plates to the bones by locking them. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-95 mm.

**Ø5.0 Locked Cannulated Screw:** Ø 5.0 mm locked cannulated screws are applied together with locked plates. It is used to fix the lower extremity plates to the bones by locking them. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-95 mm.

**Ø6.5 mm Locked Spongiosis Screw:** Ø 6.5 mm locked spongiosis screws are applied together with locked plates. It is used to fix the spongiosis structures in the distal and proximal parts of the bones of the lower extremity plates. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 30-120 mm.

**Ø6.5 Spongiosis Screw:** Ø 6.5 mm spongiosis screws are used in lower extremity anatomical plates. There are 16 mm groove and 32 mm groove options. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 30-120 mm.

**Ø4.5 mm Cortical Screw:** Ø 4,5 mm cortical screws can be unlocked in the combi-hole structures of locked plates and are applied for compression purposes. It is used to fix the lower extremity plates to the bones. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-90 mm.

**Ø3.5 mm Locked Screw:** Ø 3.5 mm locked screws are applied together with locked plates. It is used to fix upper extremity plates, pelvis plates and lower extremity plates such as distal tibia and fibula. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-70 mm.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Sti.  
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye  
**Tel:** 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

# ZİMED MEDİKAL

## LOWER-UPPER EXTREMITY BONE PLATES INSTRUCTIONS FOR USE

Submitted to the attention of Surgeons

(EN)

Document No: IFU.01.01/00

14.06.2022

**Ø3.5 Locked Cannulated Screw:** Ø 3.5 mm locked cannulated screws are applied together with locked plates. It is used to fix the upper extremity plates to the bones by locking. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-55 mm.

**Ø 4.0 mm Locked Spongiosis Screw:** Ø 4.0 mm locked cannulated spongiosis screws are applied together with locked plates. It is used to fix the spongiosis structures in the distal and proximal parts of the bones of the lower extremity plates. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 30-50 mm.

**Ø 4.0 mm Spongiose Screw:** Ø 6.5 mm spongiosis screws are used in lower extremity anatomical plates. It has a 1/3 mm groove. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 30-50 mm.

**Ø3.5 mm Cortical Screw:** Ø 3.5 mm cortical screws can be unlocked in the combi-hole structures of locked plates and are applied for compression purposes. It is used to fix upper extremity plates, pelvis plates and lower extremity plates such as distal tibia and fibula. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-120 mm.

**Ø2.4 mm cortical screw:** Ø 2.4 mm cortical screws are applied together with locked plates. It is used to apply more screws to less area of the plates such as Distal Radius, Distal Ulna, Radial Head, Distal Clavicle and to provide beter fixation. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 10-30 mm.

**Ø2.4 mm locked screw:** Ø 2.4 mm locked screws are applied together with locked plates. It is used to apply more screws to less area and fix the plates such as Distal Radius, Distal Ulna, Radial Head, Distal Clavicle. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 10-30 mm.

**Ø3.5 mm Cannulated Screw:** Ø3.5 mm locked cannulated screws are applied together with locked plates. It is used for locking the upper extremity plates to the bones. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-55 mm.

**Ø4.5 mm Cannulated Screw:** Ø4.5 mm locked cannulated screws are applied together with locked plates. It is used to fix the upper extremity plates to the bones by locking. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 12-55 mm.

**Ø 5.0 mm Locked Screw Self Drilling:** It is used to fix the lower extremity area minimally invasive plaques. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 18-115 mm.

**Ø7.3 mm Locked Cannulated Screw:** It is used to fix the locked proximal lateral femoral head in the proximal region of the plates. There are various lengths between 50-120 mm. There are 16 mm and 32 mm thread options. It is produced from Ti6Al4V ELI (ASTM F 136) material.

**Do not use the product rather than the indications stated here.**

**PATIENT POPULATION:** Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

**CONTRAINDICATIONS:** Sepsis, malignant primary or metastatic tumors, material tenderness, damaged vasculature, pertrochanteric fractures or trochanteric fractures with subtrochanteric enlargement (31-A3), overweight patients, patients with insufficient skin

tissue, bone or neurovascular condition, damaged tendon, patients with high levels of activity, patients with diabetes, patients with osteoporosis, patients with blood circulation problems, patients with material sensitivity, psychologically or physiologically incompetent patients, elderly patients, active and latent infection.

### **WARNINGS:**

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products from different manufacturers must not be used together.
4. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
5. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning and must inform the patient.
6. These products are for single use only. Do not re-use the product.
7. The products are anatomically designed. Plate bender should not be used for anatomical products.

### **SITUATIONS THAT MAY AFFECT INTENDED PERFORMANCE:**

#### **Preoperative:**

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
2. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.
3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.
4. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.
5. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.
6. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

#### **Operation**

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.
2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.
3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.
4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

#### **Postoperative**

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.
2. The patient must be warned that the implant cannot replace a normal and healthy bone, the implant can break or be damaged as a result of some movements and traumas, its expected life is limited, and replacement can be necessary in the future.
3. Abrasion of the components can result in the increase of loosening and damaging of the bone.
4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending of cracking of components.
5. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.
6. Recommendations of surgeons in the postoperative period are of utmost importance.
7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

#### **SIDE EFFECTS:**

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.
2. Damaging of blood vessels or hematomas can be seen.
3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
4. Unwanted lengthening or shortening of the extremity can occur.
5. Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss o fixation, lack of union or excessive weight
6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
8. Trochanteric non-union may occur due to inadequate re-attachment and/or premature weight bearing.
9. Trochanteric avulsion may occur as a result of excessive muscle tension, premature weight loss, or accidental intraoperative attenuation.
10. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
11. Embolism can occur in patients who have extreme sedentary lifestyle.

#### **REMOVAL OF THE IMPLANT:**

##### **Implants Put on the Market as Non-Sterile:**

Products put on the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in. all Lower-Upper Extremity Bone Plates products and surgical instruments used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

**Autoclave Sterilization:** Non-sterile products put on the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Sti.

Aydınlı Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye

Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

# ZİMED MEDİKAL

## LOWER-UPPER EXTREMITY BONE PLATES INSTRUCTIONS FOR USE

Submitted to the attention of Surgeons  
(EN)

Document No: IFU.01.01/00  
14.06.2022

AUTOCLAVE STERILIZATION	
Parameters	Pre-Vacuum Autoclave
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying time	15 minutes

**Note: Exposure time is the period to be applied starting from the moment that the autoclave reaches the exposure temperature.**

The following table shows the autoclave types as compared according to certain sterilization parameters.

TYPE OF PRODUCT	RECOMMENDED STERILIZATION METHOD
Lower-Upper Extremity Bone Plates	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vakum Otoklav): <b>If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)</b>

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

#### SHELF LIFE:

There is no shelf life for non-sterile products.

#### DISPATCH OF THE PRODUCTS:

Lower-Upper Extremity Bone Plates presented as non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

#### STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity.

**SEVERE EVENTS:** Any severe events arising in relation with the implant must be notified to the producer and used and/or authorized entity of the member state that the patient is a citizen of.

**LABEL INFORMATION:** Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

**(See: Description of Symbol)**

SYMBOL DESCRIPTION	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single-Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Used only with prescription.
	Implants are for single use
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution
	Non sterile

**Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.**

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Sti.  
Aydınlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye  
**Tel:** 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

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