www.zimed.com.tr

zimed[®]

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(EN) PRODUCT DEFINITION:

It includes plates, screws, and hand tools that will enable the application of these products to ensure the anatomical alignment of the broken bones, to keep them together, and to provide splinting with the ability to ensure the early mobility of the patient. The products are non-sterile and the implants in the set are for single use. See the product catalog for reference numbers and names of products.

DHS DCS Plate Systems are made of titanium alloy (ISO 5832-3, ASTM F 136).

DHS Plate 135°: Locked DHS Plate 135° is indicated for pertrochanteric fractures of Types 31-A1 and 31-A2, intertrochanteric fractures of type 31-A3, Basilar neck fracture 31-B, Subtrochanteric fractures and femoral neck fractures. It is used together with Ø5mm locked screw, Ø4.5 cortical screw and Ø6.5 spongiosis screw. Size options are available from 3 to 12 holes.

DCS Plate 95°: Locked DCS Plate 95° is indicated for Type 33-A fractures (extra-articular, supracondylar) Type 33-C fractures (fully articular fractures), intercondylar fractures of the femur. It is used with Ø 5mm locked screw, Ø4.5 cortical screw and Ø6.5 spongiosis screw. Size options are available from 3 to 12 holes.

DHS/DCS Lag Screw: Lag Screw is applied together with DHS and DCS plates. There are different length options between 50 mm and 120 mm.

Lag Compression Screw: Lag Compression Screw is applied together for the compression of Lag screws.

Ø6.5 mm Spongiosis Screws: Formed in a low profile structure, they can be used alone or used for fixation with plate. It can be used for plate fixation in fracture fixations, as well as for free fragments and reduction in the spongiosis region or alone for compression. With full groove and half groove options, the adherence area where compression is desired varies according to the physician's request. The head and neck of the screw are designed to fit into the hole slots in the plates. It is possible to use screw and tapped ferrule implant together and it has partial groove and full groove options and self-tapping features. There are 19 lengths with a diameter of Ø6.5 mm with a 16 mm groove and between 30-120 mm lengths, increasing by 5 mm. There are 17 lengths with a diameter of Ø6.5 mm, with a groove of 32 mm, and between 40-120 mm in length with an increase of 5 mm. There are 19 lengths with a diameter of Ø6.5 mm, fully grooved, between 30-120 mm in length, with an increase of 5 mm.

Ø4.5 mm Cortical Screws: They are used in the cortical region as well as in plate fixation in fracture fixations. It provides full harmony with the shape on the head and the dynamic hole in the plate, and in this way, the screw can be completely in contact with the plate and can be positioned freely. Thanks to its self-tapping feature, it allows fast and easy application. With a diameter of Ø4.5 mm, it increases by 2 mm between 12-60 mm lengths and has a total of 31 length options with an increase of 5 mm between 60-90 mm lengths.

Ø5.0 mm Low Profile Locking Screws: They are used in the cortical region as well as in plate fixation in fracture fixations. Thanks to the grooves on the head, it provides a fully compatible locking with the locking plate and thus completely eliminates the risk of screw loosening. The head and neck parts are designed to fit perfectly into the hole slots on the plates. There are 29 length options with a diameter of Ø5.0 mm, with a height increase of 2 mm between 12-50 mm and a height increase of 5 mm between 50-95 mm.

Do not use the product rather than the indications stated here.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

ZİMED MEDİKAL

DHS&DCS PLATE SYSTEM INSTRUCTIONS FOR USE

Submitted to the attention of Surgeons

(EN)

Document No: IFU.01.02/00

14.06.2022

4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.
2. The patient must be warned that the implant cannot replace a normal and healthy bone, the implant can break or be damaged as a result of some movements and traumas, its expected life is limited, and replacement can be necessary in the future.
3. Abrasion of the components can result in the increase of loosening and damaging of the bone.
4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending or cracking of components.
5. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.
6. Recommendations of surgeons in the postoperative period are of utmost importance.
7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

SIDE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.
2. Damaging of blood vessels or hematomas can be seen.
3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
4. Unwanted lengthening or shortening of the extremity can occur.
5. Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss of fixation, lack of union or excessive weight
6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
8. Trochanteric non-union may occur due to inadequate re-attachment and/or premature weight bearing.
9. Trochanteric avulsion may occur as a result of excessive muscle tension, premature weight loss, or accidental intraoperative attenuation.
10. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
11. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT:

Implants Put on the Market as non-Sterile:

Products put on the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in. all DHS&DCS Plate System products and surgical hand tools used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 44 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

TR
ENGES
IT

FR

RU

ZİMED MEDİKAL
DHS&DCS PLATE SYSTEM INSTRUCTIONS FOR USE
Submitted to the attention of Surgeons
(EN)

Document No: IFU.01.02/00
14.06.2022

Autoclave Sterilization: Non-sterile products put on the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

AUTOCLAVE STERILIZATION	
Parameters	Pre Vacuum Autoclave
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying time	15 minutes

Note: Exposure time is the period to be applied starting from the moment that the autoclave reaches the exposure temperature.

The following table shows the autoclave types as compared according to certain sterilization parameters.

TYPE OF PRODUCT	RECOMMENDED STERILIZATION METHOD
DHS&DCS Plate System	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vakum Otoklav): <i>If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)</i>

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

SHELF LIFE:
There is no shelf life for non-sterile products.

DISPATCH OF THE PRODUCTS:
DHS&DCS Plate System presented as non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:
Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the producer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.
(See: Description of Symbol)

SYMBOL DESCRIPTIONS	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single-Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Used only with prescription.
	Implants are for single use
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution

	Non sterile
--	-------------

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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