[www.zimed.com.tr](http://www.zimed.com.tr)

**zimed®**

**ELASTİK ÇİVİLER  
KULLANIM KİLAVUZU  
INSTRUCTIONS FOR USE  
OF ELASTIC NAIL**

**TR**

**ENG**

**ES**

**IT**

**FR**

**RU**

*innovation for move*

# ZİMED MEDİKAL

## ELASTİK ÇİVİLER KULLANIM KİLAVUZU

### Uzman Hekimin İlgisine (TR)

Doküman No: IFU.04.02/00  
14.06.2022

#### (TR) ÜRÜN TANIMI:

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**(EN) PRODUCT DEFINITION:**

Titanium Elastic Nail is designed to be used in pediatric patients for the fixation of diaphyseal fractures where the medullar canal is narrow or where flexibility of the implant is important. At the same time, it is also indicated for the treatment of epiphyseal and metaphyseal fractures such as radial neck fractures and fixation of long and thin bones such as carpal or tarsal bones. Elastic nail is 440 mm long, and diameter options include Ø 2, (2,5), 3, (3,5), 4, (4,5), 5 and (5,5) mm. Different colors represent different diameter measures including yellow, blue, pink, green and brown. Elastic Nail is used with an end cap suiting the diameter measures including (2, (2,5), 3, (3,5), 4, (4,5), 5 and (5,5) mm). Titanium Elastic Nail is made from Ti6Al7Nb.

**INTENDED USE AND INDICATIONS:**

Titanium Elastic Nail is designed to be used in pediatric patients for the fixation of diaphyseal fractures where the medullar canal is narrow or where flexibility of the implant is important. At the same time, it is also indicated for the treatment of epiphyseal and metaphyseal fractures such as radial neck fractures and fixation of long and thin bones such as carpal or tarsal bones.

**Please do not use this product other than the indications stated here.**

**PATIENT POPULATION:** Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

**CONTRA-INDICATIONS:**

Active infections.  
Poor bone quality or amount, disorders in the medullar canal, poor blood circulation, and previous infections will delay healing  
Manifest or suspected occult infection or manifest local inflammation in or around the affected site  
Dangers related to vascularization that could hinder adequate supply of blood to the fracture or operation site,  
Existing fracture fixation device implanted before  
Bone deformation that had appeared in the past,  
Hypovolemia, hypothermia, coagulopathy,  
Mental disorders preventing cooperation during the rehabilitation regime,  
Obligation to pass through the active epiphyseal plaques,  
Low sub-trochanteric fractures,  
Femoral shaft fractures,  
Isolated or combined medial femoral neck fractures,  
Poor bone stock will have negative effects on the disease, can lead to infections and change the priorities.

**WARNINGS:**

- Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
- In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
- Products of different manufacturers must not be used together.
- Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
- Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning and must inform the patient.
- These products are for single use only. Do not re-use the product.

# ZİMED MEDİKAL

## INSTRUCTIONS FOR USE FOR ELASTIC NAILS

### Submitted to the Attention of Surgeons (EN)

Document No: IFU.04.02/00

14.06.2022

7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

**ADVERSE EFFECTS:**

- Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.
- Damaging of blood vessels or hematomas can be seen.
- Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
- Unwanted lengthening or shortening of the extremity can occur.
- Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss of fixation, lack of union or excessive weight.
- As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
- Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
- Trochanteric nonunion can happen with the reason of inadequate re-union and/or load-bearing in the early period.
- Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
- Embolism can occur in patients who have extreme sedentary lifestyle.

**REMOVAL OF THE IMPLANT:**

The period of the implant within the patient's body is 12 to 18 months and must be removed according to the doctor's decision.

**STERILIZATION:**

**Implants Put in the Market as Non-Sterile:** Products put in the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in.

Elastic Nail System products and surgical tools used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

**Autoclave Sterilization:** Non-sterile products put in the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

AUTOCLAVE STERILIZATION	
Parameters	Autoclave with pre-vacuum
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying Time	15 minutes

**Note: Exposure time means the period to be applied starting from the moment that the autoclave temperature reaches the exposure temperature.**

The table below shows autoclave types for specific sterilization parameters.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Şti.

Aydınlar Mah. 03070 Cad. No:4 Sehitkâmil 27580 Gaziantep-Turkey

Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

**ZİMED MEDİKAL**  
**INSTRUCTIONS FOR USE FOR ELASTIC NAILS**  
Submitted to the Attention of Surgeons  
(EN)

Document No: IFU.04.02/00  
14.06.2022

PRODUCT TYPE	RECOMMENDED STERILIZATION METHOD
Elastic Nail System	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vacuum Autoclave): <b><i>If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)</i></b>

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

**SHELF LIFE:**  
Non-sterile products do not have shelf lives.

**DISPATCH OF THE PRODUCTS:**  
The Elastic Nail System presented as non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrinked together with instructions for use and label.

**STORAGE AND MAINTENANCE CONDITIONS:**  
Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity.

**SEVERE EVENTS:** Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

**LABEL INFORMATION:** Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.  
(See: Description of Symbol)

SYMBOL DESCRIPTIONS	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single Device Identification Number)
	Electronic Instructions for Use
	Manufacturer Details
	Production Date
	Used only with prescription
	Implants are for single use only.
	Do not use the damaged package.
	Keep away from direct sunlight.
	Keep away from precipitation.
	Fragile object: Handle with care
	Temperature limits
	Caution
	Non sterile

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

**MANUFACTURER:** Zimed Medikal San. ve Tic. Ltd. Şti.  
Aydınlar Mah. 03070 Cad. No:4 Sehitkamil 27580 Gaziantep-Turkey  
Tel: 0 342 238 44 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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