www.zimed.com.tr

zimed[®]

**KALÇA EKLEM PROTEZLERİ
KULLANIM KILAVUZU**
*INSTRUCTIONS FOR USE
OF HIP JOINT PROSTHESES*

TR

ENG

ES

IT

FR

RU

innovation to move

TR

ENG

ES

IT

FR

RU

**Operasyon**

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CE 1984

ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF HIP JOINT PROSTHESES

Submitted to the Attention of Surgeons

(EN)

Document No: IFU.07/00
14.06.2022

(EN) 1. PRODUCT GROUP: Hip Joint Prostheses

Uncemented Acetabular Cup acts as the acetabular cup in patients with tumors and abraded acetabular area. It is made from Ti 6Al4V ASTM F136 (ISO 5832-3) raw material in anatomic structure, and has no right/left differentiation. Applications with HA, Dual and Plasma Porous coatings are manufactured. Uncemented Acetabular Cup has 11 sizes varying diameters between 42-62 mm. Holes of the Uncemented Acetabular Cup are designed so as to have 2 holes in diameters between 42 and 50 mm and 3 holes in diameters between 52 and 62 mm. Uncemented Acetabular Cup is sterilized with gamma irradiation and EO sterilization method.

Acetabular Liner is made from UHMWPE raw material complying with the ISO 5834-2 standard. Acetabular Liner becomes functional as a component of the joint after being mounted with the Acetabular Cup and the Femoral Head, one inside another. Acetabular liner has 11 sizes varying diameters between 42-62 mm. It is used with Heads with diameters Ø22, Ø28, Ø32 and Ø36 mm. Use of the acetabular liner and the femoral head is given in the table below.

		Acetabular Liner										
		42	44	46	48	50	52	54	56	58	60	62
Femoral Head	22											
	28											
	32											
	36											

Acetabular Screw is made from Ti6Al4V ELI raw material complying with ASTM F136 (ISO 5832-3) standard. Acetabular Screw provides for the fixation of Acetabular Cup in acetabular area.

Cemented Acetabular Cup is made from UHMWPE raw material complying with ISO 5834-2 standard. It has no right/left differentiation. Cemented acetabular cup is used together with heads with diameters of Ø28 mm r Ø36 mm head.

		Cemented Acetabular Cup											
		Ø	42	44	46	48	50	52	54	56	58	60	62
Femoral Head	28												
	36												

Femoral Head is produced with the diameters of Ø22, Ø28, Ø32 and Ø36 mm, and is made from both CoCrMo Alloy (ISO 5832-4) and Implant Steel (ISO 5832-1) raw material and has anatomic structure, Femoral Head provides for the fixation of and connection between the Femoral Stem and Acetabulum.

Bipolar Head: The outer portion is made from both CoCrMo Alloy (ISO 5832-4) and Implant Steel raw material (ISO 5832-1); while the inner portion is made from UHMWPE raw material. Bipolar Head products are used together with Heads with diameters of 22 and 28 mm. Ways of using are given in the table below.

		Bipolar Head											
			38	40	42	44	46	48	50	52	54	56	58
Femoral Head	22												
	28												

Unipolar Head is made from Implant Steel (ISO 5832-1) raw material. It has no right/left differentiation. It has several diameter sizes in the range of 40 mm – 56 mm.

Stem with Calcar Support is made from CoCrMo Alloy (ISO 5832-4) and Implant Steel (ISO 5832-1) raw material. HA, Dual and Porous coating is used in uncemented applications. It substitutes for bone in fractures and tumors of the femoral head area and provides support for the body.

Femoral Stem is made from CoCrMo Alloy (ISO 5832-4) and Implant Steel (ISO 5832-1) raw materials for cemented applications and from Ti6Al4V ELI in accordance with STM F136 (ISO 5832-3) standard for uncemented (with HA, Dual and Porous coating) applications to substitute for bone in fractures and tumors in femoral head and neck areas and to provide support to the body.

Straight Stem is designed for use in the elderly with Collum Femoris fractures and in patients who needs total hip prosthesis. It is made from implant Steel raw material (ISO 5832-1), and is applied only with cement.

Cement Stopper is made from UHMWPE raw material complying with the ISO 5834-2 standard. It prevents the flow of the cement within the bone and helps for holding of the prosthesis.

Distal Centralizer is made from UHMWPE raw material. Distal Centralizer ensures centralization of the prosthesis in the bone, and provides for accurate positioning of the prosthesis. Stems fitting the measures of the centralizer are defined in the following table.

Distal Centralizer		Cemented Femoral Stem (Size)	Stem supported with Cemented Calcar (Size)
	8,5 mm	7 mm, 8 mm	9 mm, 11 mm
	11,5 mm	9 mm, 10 mm	13 mm
	13,5 mm	11 mm, 12 mm	

INTENDED USE AND INDICATIONS: Hip Joint Prostheses

It is indicated in degenerative osteoarthritis, displaced intracapsular femoral neck fractures, treatment inter-trochanteric fractures, in the elderly patients with collum femoris fractures, traumatic arthritis, slipped femoral head epiphysis, united hip and dystrophic changes, primarily in femur with poor trochanter minor and calcar support, inflammatory diseases such as rheumatoid arthritis causing intra-articular damage, hip calcification arising from congenital hip dislocation and poor development of the hip joint, avascular necrosis (softening of bone in relation with the functional loss of vessels supplying the head part of the femur and consequential destruction of the hip joint), previous failed surgical interventions or failure in obtaining satisfactory stabilization in flexion and definitive diagnosis of any one or more conditions mentioned above, and in primary or revision surgeries in cases where other treatments or devices have proven unsatisfactory. At the same time, in proximal femoral fractures involving the femoral head and neck that correction with other techniques is not possible, in endoprotheses, in femur osteotomy or Girdlestone resection, slipped hip after fractures and in correction of deformations.

Please do not use the products rather than the indications stated here.

CONTRA-INDICATIONS:

Conditions that block/tend to block adequate support of the implant or prevent use of an implant with suitable measures, including, for example:

- o Limited blood supply;
- o Poor bone quality and amount; for example, osteoporosis or other metabolic disorders that weaken bone formation and osteomalacia;
- o Infections or other conditions leading to advanced bone resorption.

Mental or nervous disorders deteriorating the ability or and motivation of the patient to limit their physical activities.

Patients with active or suspicious infection in the area surrounding the hip joint, Physical conditions and movements resulting in the overloading of implants, for example, Charcot joints, missing muscles, multiple joint mobility disorders, etc.

Skeletal immaturity,

Patients with any vascular insufficiency, muscular atrophy or serious neuromuscular disorders,

Patients with known medium- or serious renal insufficiency,

Patients with suppressed immune systems because of AIDS or use of corticosteroid in high dosages,

Physical disorders such as morbid obesity or limited motility, Failure of the patient in following the instructions for postoperative care, Contra-indications can be relative or definitive. Patients must be examined very carefully, and alternative procedures must be used if possible, including for example, nonsurgical treatment, osteosynthesis, arthrodesis, femoral osteotomy, pelvic osteotomy, resection arthroplasty, half arthroplasty and others.

Conditions involving high risks include the following: osteoporosis, metabolic disorders preventing bone formation and osteomalacia.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the hip joint prosthesis and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if s/he had previous operations or not etc. into consideration.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.

2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.

3. Products of different manufacturers must not be used together.

4. Components that can be used together according to EN ISO 21534 Annex-C:

CoCrMo Alloy (ISO 5832-4)/ Titanium Alloy (ISO 5832-3)

CoCrMo Alloy (ISO 5832-4)/ CoCrMo Alloy (ISO 5832-4)

Implant Steel (ISO 5832-1)/Titanium Alloy (ISO 5832-3)

Implant Steel (ISO 5832-1)/Implant Steel (ISO 5832-1)

Components that are incompatible with each other:

Implant Steel (ISO 5832-1)/ CoCrMo Alloy (ISO 5832-4)

5. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.

6. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning, and must inform the patient.

7. These products are for single use only. Do not re-use the product.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.

2. Selection of the prosthesis suiting the fracture type is of utmost importance. The surgeon must select the prosthesis suiting the image in the x-ray before the operation. Use of uncemented femoral stem is recommended for 31A1 fracture type.

3. The surgeon must warn the patient that the prosthesis cannot replace a normal healthy bone, and the prosthesis can break or be damaged as a result of some movements or traumas, the expected lifespan is limited, and replacement in future can be required.

4. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.

5. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlık Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

ZİMED MEDİKAL

INTRUCTIONS FOR USE OF HIP JOINT PROSTHESES

Submitted to the Attention of Surgeons
(EN)

Document No: IFU.07/00
14.06.2022

6. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.

7. Incorrect selection or positioning of components can result in insufficient femoral range of motion because of femoral squeeze and peri-articular calcification.

8. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

9. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.

2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.

3. Correct selection of the shell and stem based on the length of the neck is very important. Muscular laxity and/or incorrect placement of components can lead to loosening, partial displacement and/or break of implants. Longer neck length and placement with inward inclination will increase the tension that should be borne by the stem. Components must be tightly implanted in place using surgical hand tools and must be checked with fluoroscopy in case of suspicion.

4. As it is not possible to see the internal tension of a prosthesis implanted before, bending or fractures in the early period can be seen; therefore, prosthesis used before must never be re-used.

5. During the implantation, femoral or acetabular perforation or fractures can occur because of trauma or overloading, particularly in the presence of poor bone structure.

6. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. The patient must be warned that the implant cannot replace a normal and healthy bone, the implant can break or be damaged as a result of some movements and traumas, its expected life is limited, and replacement can be necessary in the future.

3. Abrasion of the components can result in the increase of loosening and damaging of the bone.

4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending of cracking of components.

5. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.

6. Recommendations of surgeons in the postoperative period are of utmost importance.

7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarct can occur.

2. Damaging of blood vessels or hematomas can be seen.

3. Delayed wound healing: Deep wound infections (early or late) that might require the removal of the prosthesis can occur. Arthrodesis of the relevant joint or amputation of the extremity can be required in rare cases.

4. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.

5. Unwanted lengthening or shortening of the extremity can occur.

6. Fatigue fractures of the components of prosthesis can arise from trauma, exhausting physical activity, misalignment, defective implantation, service period, loss of fixation, nonunion or overloading.

7. Trochanteric nonunion can occur because of insufficient re-connection or early load-bearing, and collapse of the prosthesis can occur.

8. Trochanteric avulsion can occur as a result of excessive muscular tension, early weight loss or inadvertent in intra-operative weakening.

9. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.

10. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.

11. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.

12. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body can change according to the daily activity of the patient. Product life of this implant has been determined as 15 years. Implants can be removed according to the surgeon's decision.

STERILIZATION:

Implants are submitted as sterile. Sterilization method will be shown on the packaging. Minimum 25 kGy gamma irradiation is applied to all the implants sterilized with gamma irradiation. Sterile products will remain so unless the integrity of the package is broken. The products are presented as sterile. **Re-sterilization** of the product is not recommended.

Please inform the manufacturer about any products with impaired sterility although offered as sterile.

SHELF LIFE:

The shelf life for **Hip Joint Prostheses** sterilized with Gamma sterilization is 5 years, while the same for those sterilized with ethylene oxide is 3 years.

Please check the expiry date. Please contact the manufacturer for products expired.

Non-sterile products do not have shelf lives.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference

to the ISO 639-1 standard in the languages of the Community.
(See: **Description of Symbol**)

DESCRIPTION OF SYMBOLS	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single-Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation.
	Sterilized with ethylene oxide
	Used only with prescription.
	Implants are for single use
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlık Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF HIP JOINT PROSTHESES

Submitted to the Attention of Surgeons
(EN)

Document No: IFU.07/00
14.06.2022

	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

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		Revestimiento Acetabular										
		42	44	46	48	50	52	54	56	58	60	62
Cabezal Femoral	22	✓	✓	✓								
	28				✓	✓	✓	✓	✓	✓	✓	✓
	32						✓	✓	✓	✓	✓	✓
	36								✓	✓	✓	✓

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		Copa Acetabular Cementada											
		Ø	42	44	46	48	50	52	54	56	58	60	62
Cabezal Femoral	28	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	36					✓	✓	✓	✓	✓	✓	✓	✓

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