www.zimed.com.tr

zimed[®]

RADYOLUSENT FİKSATÖRLER
KULLANIM KILAVUZU
INSTRUCTIONS FOR USE
OF RADIOLUCENT FIXATORS

TR

ENG

ES

IT

FR

RU

innovation to move

ZİMED MEDİKAL

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ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF RADIOLUCENT FIXATORS

Submitted to the Attention of Surgeons

(EN)

Document No: IFU.05.03/00
14.06.2022

(EN) PRODUCT DEFINITION:

Radiolucent Wrist Fixator: The wrist fixator with dynamic-angle clamp is used in the partial osteotomy of distal radius deformity, and acute correction by filling the space with graft of gradual correction of the space with hemicallosis osteotomy. It can be used in distal radial and/or ulnar fractures and unstable intra-articular distal radial fractures (by fixing the fracture fragments with other implants). The ball-joint structure allows angulated locking.

INDICATIONS:

Radiolucent fixator is indicated in fixation of Colles fractures and distal radius osteotomy procedures. It has not been designed for use in the spine. Radiolucent fixator is a radiolucent device designed to help in the repair of unstable distal radius fractures. This structure is designed to help the surgeon in obtaining the initial reduction while the fixator is in place and then to adjust the planes required only for correction.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

CONTRA-INDICATIONS:

Active infections, manifest or occult infection or manifest local inflammation in or around the affected site, existing fracture fixation device implanted before, existing bone deformation, hypovolemia, hypothermia, coagulopathy, obesity, advanced osteoporosis, poor bone quality or amount, disorders in the medullar canal, and poor blood circulation are among the contra-indications. It is contra-indicated in patients with mental or physiologic disorders who are unwilling or unable to comply with postoperative care instructions. It is not recommended in patients with severe osteoporosis, who are HIV-positive or with immunosuppression, or poorly-controlled diabetes mellitus. Tests must be carried out before placement of the implant in patients with foreign body hypersensitivity or suspected sensitivity against materials.

WARNINGS:

- Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
- In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
- Products of different manufacturers must not be used together.
- Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
- Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning and must inform the patient.
- These products are for single use only. Do not re-use the product.
- Products must never be rinsed with ultrasonic washing with detergent. Ultrasonic washing with detergent can result in the removal of the eloxal coating of external fixator products.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

- The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
- Surgeon must also inform the patient about any other risks s/he can deem required to be acknowledged to the patient.

3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, moment of resistance and/or erosion characteristics of the implant. These can transform into invisible internal tension and can result in breaking of the product.

4. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.

5. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

6. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative:

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.

2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.

3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.

4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. Long-term periodic follow-up is recommended to observe the positions and conditions of the prosthetic components and the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending or cracking of components.

3. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.

4. Recommendations of surgeons in the postoperative period are of utmost importance.

5. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.

2. Damaging of blood vessels or hematomas can be seen.

3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.

4. Unwanted lengthening or shortening of the extremity can occur.

5. Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss of fixation, lack of union or excessive weight

6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.

7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.

8. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.

9. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body is minimum 6 weeks and can be extended to 1 year or more in difficult situations; afterwards, it must be removed based on the surgeon's decision.

STERILIZATION:

Radiolucent fixators are put in the market as sterile. Sterile products will remain so unless the integrity of the packaging is broken. Two types of sterilization methods, gamma irradiation and ethylene oxide, are applied. Minimum 25 kGy gamma irradiation dosage is applied to all the implants sterilized with gamma irradiation.

Please notify the manufacturer about the product with impaired sterility if put into the market as sterile.

SHELF LIFE:

Shelf life for Radiolucent fixators sterilized with gamma irradiation is 5 years; while the same for those sterilized using ethylene oxide is 3 years.

Please check the expiry date. Contact the manufacturer for products expired.

DISPATCH OF THE PRODUCTS:

Radiolucent Fixators presented as sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the country code in ISO 639-1 standard in the languages of the Community. **(See: Description of Symbols)**

DESCRIPTION OF SYMBOLS	
	Logo of the Firm
	"1984" Notified Body Number
	Medical Device

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlık Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web:www.zimed.com.tr E-mail: info@zimed.com.tr

ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF RADIOLUCENT FIXATORS

Submitted to the Attention of Surgeons
(EN)

Document No: IFU.05.03/00
14.06.2022

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

	Lot Number
	Product Reference/Catalog Number
	UDI Number (Single Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation
	Sterilized with ethylene oxide
	Used with prescription
	Implants are for single use only.
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

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