www.zimed.com.tr

zimed®

**REVİZYON DİZ PROTEZLERİ
KULLANIM KİLAVUZU
INSTRUCTIONS FOR USE
OF REVISION KNEE PROSTHESES**

TR

ENG

ES

IT

FR

RU

innovation for move

TR

ENG

ES

IT

FR

RU

ZİMED MEDİKAL

REVİZYON DİZ PROTEZLERİ KULLANIM KİLAVUZU

Uzman Hekimin İlgisine (TR)

Doküman No: IFU.09/00
14.06.2022

(TR)

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(EN)

1. PRODUCT GROUP: Revision Knee Prostheses

Tibial Plateau (Tibial Component): It is used in the revision of tibial defects. It is made from CoCrMo Alloy complying with ISO 5832-4. It has 6 different sizes including 2, 3, 4, 5, 6 and 7 mm.

Tibial Insert: It is used as a connector in revision knee prostheses with the purpose of providing wider flexion range in patients with eligible anatomy. It is designed so as to allow full flexion movement between 140° and 150°. This internal component is made from Ultra High Molecular Weight Polyethylene (UHMWPE) complying with ISO 5834-2 standard with high resistance. Thickness options include 9, 11, 13, 15, 17, 19, and 21 mm to help the establishment of soft tissue balance in the knee joint.

Flat Stem: It is made from Ti6Al4V Eli raw material complying with ISO 5832-3 standard. Flat stems are used to reduce the flexion force in the knee in revision procedure of knee prostheses and to distribute this force onto the femoral and tibial axes. Diameter options include 9, 10, 11, 12, 13, 14, 15 and 16 mm, and proximal length options include 60, 100, and 130 mm.

Offset Stem: It is made from Ti6Al4V Eli material complying with ISO 5832-3 standard. Offset stems are used with the purpose of positioning the femoral and tibial components in the best way to comply with the patient's anatomy. Diameter options include 9, 10, 11, 12, 13, 14, 15 and 16 mm, and proximal length options include 60, 100, and 130 mm.

Femoral Component: It is used in the revision of the femoral defect. It is made from CoCrMo alloy complying with ISO 5832-4 standard. It has anatomical structure, and it differs for right and left sides. It has 6 different sizes as 2, 3, 4, 5, 6 and 7 mm.

Distal Femoral Wedge: It is made from Ti6Al4V material complying with ISO 5832-3 standard. It has anatomical structure and does not differ for right and left sides. It is used to fill up the femoral knee defects. Optional thicknesses are 5 and 10. It is fixed onto the components with the help of screws.

Posterior Condylar Wedge: It is made from Ti6Al4V material complying with ISO 5832-3 standard. It has anatomical structure and does not differ for right and left sides. It is used to fill up the femoral knee defects. Optional thicknesses are 5 and 10. It is fixed onto the components with the help of screws.

Tibial Half Wedge: It is made from Ti6Al4V Eli material complying with ISO 5832-3 standard. It has anatomical structure and does not differ for right and left sides. It is used to fill up the tibial knee defects. Optional thicknesses are 5 and 10. It is fixed onto the components with the help of screws.

Tibia Insert Screw: It is used to fix the tibial insert. It is made from Ti6Al4V Eli material complying with ISO 5832-3 standard. It has one size option of Ø 5x33 mm.

Wedge Screw: It is used to fix the wedge onto the femoral component. It is made from Ti6Al4V Eli material complying with ISO 5832-3 standard. It has two types as the tibial and femoral wedge screws. It has two size options of 5 and 6 mm at the diameter of Ø 5 mm.

INTENDED USE AND INDICATIONS:**Revision Knee Prostheses:**

The need for revision in primary total knee joint prosthesis procedures arises in time because of loosening of one or more components because of friction, from wear and tear, and infections; poor bone growth on uncemented implant surface, or patellar problems. Risk of infection in total knee prosthesis is less than 1%; however, knee revision operation can be required if infection occurs. Revision Knee Prosthesis is used in cases with severe knee pain and in disabilities with various reasons. Conditions leading damage or disability in knee include the following:

Rheumatoid arthritis, osteoarthritis, traumatic arthritis, and polyarthritis,
Collagen disorders and/or avascular necrosis of the femoral condyle,
Valgus, varus or flexion deformities of medium level,
Previously failed surgical interventions or failure in obtaining satisfactory stability in flexion.

Loosening of knee prosthesis because of traumatic fractures or displacement or failure in obtaining stability of the knee.

Do not use the product rather than the indications stated here.

ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF REVISION KNEE PROSTHESES

Submitted to the Attention of Surgeons (EN)

Document No: IFU.09/00

14.06.2022

CONTRA-INDICATIONS:

Contra-indications of Revision Knee Prosthesis and undesirable conditions that might occur:
Local/systemic infections involving the joint or past infections in the history that might affect the prosthesis,

Poor bone stock in femoral or tibial surfaces,
Skeletal immaturity,
Neuropathic arthropathy,
Osteoporosis, any muscular loss or neuromuscular diseases,
Severe instability secondary to the absence of collateral ligament integrity,
Mental disorders that might deteriorate patient compliance.

Total knee arthroplasty is contra-indicated in patient with rheumatoid arthritis (RA) or skin ulcers or repetitive destruction of the skin, because postoperative infection risk will be higher in such cases. Infection risk may be even higher in rheumatoid arthritis patients who use steroids. Time for late infections in RA patients has been reported as postoperative 24+ months.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the knee joint prostheses and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if s/he had previous operations or not etc. into consideration.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.

2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.

3. Products of different manufacturers must not be used together.

4. Components that can be used together according to EN ISO 21534 Annex-C:

CoCrMo Alloy (ISO 5832-4)/ Titanium Alloy (ISO 5832-3)
CoCrMo Alloy (ISO 5832-4)/ CoCrMo Alloy (ISO 5832-4)
Implant Steel (ISO 5832-1)/Titanium Alloy (ISO 5832-3)
Implant Steel (ISO 5832-1)/Implant Steel (ISO 5832-1)

Components that are incompatible with each other:

Implant Steel (ISO 5832-1)/ CoCrMo Alloy (ISO 5832-4)

5. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.

6. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning and must inform the patient.

7. These products are for single use only. Do not re-use the product.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:**Preoperative:**

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.

2. The surgeon must warn the patient that the prosthesis cannot replace the normal healthy bone, and the prosthesis can break or be damaged as a result of some movements or traumas, the expected lifespan is limited, and replacement in future can be required.

3. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.

4. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension and can result in breaking of the product.

5. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.

6. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

7. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative:

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.

2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.

3. As it is not possible to see the internal tension of a prosthesis implanted before, bending or fractures in the early period can be seen; therefore, prosthesis used before must never be re-used.

4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. Abrasion of the components can result in the increase of loosening and damaging of the bone.

3. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending or cracking of components.

4. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.

5. Recommendations of surgeons in the postoperative period are of utmost importance.

6. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarct can occur.

2. Damaging of blood vessels or hematomas can be seen.

3. Delayed wound healing: Deep wound infections (early or late) that might require the removal of the prosthesis can occur. Arthrodesis of the relevant joint or amputation of the extremity can be required in rare cases.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

TR

ENG

ES

IT

FR

RU

ZİMED MEDİKAL
INSTRUCTIONS FOR USE OF REVISION KNEE PROSTHESES
Submitted to the Attention of Surgeons
(EN)

Document No: IFU.09/00
14.06.2022

4. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.

5. Unwanted shortening or lengthening of the extremity can occur.

6. Fatigue fractures of the component of the prosthesis can arise from trauma, exhausting physical activity, misalignment, defective implantation, service period, loss of fixation, nonunion or overloading.

7. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.

8. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.

9. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.

10. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body can change according to the daily activity of the patient. Product life of this implant has been determined as 15 years. Implants can be removed according to the surgeon's decision.

STERILIZATION:

Implants are submitted as sterile. Sterilization method will be shown on the packaging. Minimum 25 kGy gamma irradiation is applied to all the implants sterilized with gamma irradiation. Sterile products will remain so unless the integrity of the package is broken. The products are presented as sterile. **Re-sterilization** of the product is not recommended.

Please inform the manufacturer about any products with impaired sterility although offered as sterile.

SHELF LIFE:

The shelf life for **Total Knee Prostheses** sterilized with Gamma sterilization is 5 years, while the same for those sterilized with ethylene oxide is 3 years.

Please check the expiry date. Please contact the manufacturer for products expired.

Non-sterile products do not have shelf lives.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: Description of Symbol)

DESCRIPTION OF SYMBOLS	
	Logo of the Firm
	"1984" Notified Body Number
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI Number (Single Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with Irradiation
	Sterilized with ethylene oxide
	Use with prescription
	Implants are for single use
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight

	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 44 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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