www.zimed.com.tr

zimed®

**REVİZYON KALÇA PROTEZİ
KULLANIM KİLAVUZU**
**INSTRUCTIONS FOR USE OF
REVISION HIP PROSTHESIS**

TR

ENG


ES

IT

FR

RU

innovation to move

ZİMED MEDİKAL

REVİZYON KALÇA PROTEZİ KULLANIM KİLAVUZU

Uzman Hekimin İlgisine (TR)

Doküman No: IFU.10/00
14.06.2022

(TR)

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ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF REVISION HIP PROSTHESIS

Submitted to the Attention of Surgeons (EN)

Document No: IFU.10/00
14.06.2022

(EN)

1. PRODUCT GROUP: Revision Hip Prosthesis

Monoblock Revision Stem: It has been designed to bridge the proximal bone defect in order to provide rotational fixation and rotational stability. A conic implant body and eight longwise canals provide the rotational stability and intra-operative anteversion adjustment. It is made from titanium alloy in accordance with ASTM F 136. Monoblock revision stem products have different diameters as 10,11,12,13,14,15 and 16, and one length option of 185 mm.

Modular Proximal Body: It is designed to be used together with modular stems. It is made from titanium manufactured in accordance with ASTM F136. Diameters include Ø 17, 20, and 23 and it has 5 different length options as 70, 80, 90, 100, and 110 mm.

Modular Flat Stem: It replaces the bone in cases of fractures and tumors in the neck area and provides support to the body. It has anatomic structure, and is made from titanium manufactured in accordance with ASTM F136. Diameters include Ø 12, 13, 14, 15, 16, 17, 18, 19, 20, 21, 22, 23 mm and the length is 120 mm. It is suitable for dual coating.

Proximal Locking Screw: Modular Proximal Body and Modular Revision Stem are designed to be interlocked with each other. It is made from titanium manufactured in accordance with ASTM F136. It has 5 different length options including 70, 80, 90, 100, and 110 mm.

Acetabular Cage: It is used in Hip Prostheses where the acetabular wall is weak or in acetabular revisions with bone defects. It is made from titanium alloy manufactured in accordance with ISO 5832-2, and it has right- and left-side options. Acetabular Cage has several size options in the range between 56 mm and 70 mm.

INTENDED USE AND INDICATIONS:

Revision Hip Prostheses

It is indicated in degenerative osteoarthritis, displaced intracapsular femoral neck fractures, treatment of inter-trochanteric fractures, in the elderly patients with collum femoris fractures, non-traumatic or non-induced wear of joints, **traumatic arthritis**, slipped femoral head epiphysis, united hip and dystrophic changes, primarily in femur with poor trochanter minor and calcar support, inflammatory diseases such as rheumatoid arthritis causing intra-articular damage, hip calcification arising from congenital hip dislocation and poor development of the hip joint, avascular necrosis (softening of bone in relation with the functional loss of vessels supplying the head part of the femur and consequential destruction of the hip joint), previous failed surgical interventions or failure in obtaining satisfactory stabilization in flexion and definitive diagnosis of any one or more conditions mentioned above, and in primary or revision surgeries in cases where other treatments or devices have proven unsatisfactory. At the same time, in proximal femoral fractures involving the femoral head and neck that correction with other techniques is not possible, in endoprostheses, in femur osteotomy or Girdlestone resection, slipped hip after fractures and in correction of deformations.

Please do not use the products rather than the indications stated here.

CONTRA-INDICATIONS:

Conditions that block/tend to block adequate support of the implant or prevent use of an implant with suitable measures, including, for example:

- Limited blood supply;
- Poor bone quality and amount; for example, osteoporosis or other metabolic disorders that weaken bone formation and osteomalacia;
- Infections or other conditions leading to advanced bone resorption.

Mental or nervous disorders deteriorating the ability or and motivation of the patient to limit their physical activities.

Patients with active or suspicious infection in the area surrounding the hip joint, Physical conditions and movements resulting in the overloading of implants, for example, Charcot joints, missing muscles, multiple joint motility losses, etc.

Skeletal immaturity, Patients with any vascular insufficiency, muscular atrophy or serious neuromuscular disorders,

Patients with known medium- or serious renal insufficiency,

Patients with suppressed immune systems because of AIDS or use of corticosteroid in high dosages,

Physical disorders such as morbid obesity or limited motility, Failure of the patient in following the instructions for postoperative care, Contra-indications can be relative or definitive. Patients must be examined very carefully, and alternative procedures must be used where possible, including for example, nonsurgical treatment, osteosynthesis, arthrodesis, femoral osteotomy, pelvic osteotomy, resection arthroplasty, half arthroplasty and others. Conditions involving high risks include the following: osteoporosis, metabolic disorders preventing bone formation and osteomalacia.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the hip joint prostheses and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if s/he had previous operations or not etc. into consideration.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products of different manufacturers must not be used together.
4. Components that can be used together according to EN ISO 21534 Annex-C:
 - CoCrMo Alloy (ISO 5832-4)/ Titanium Alloy (ISO 5832-3)
 - CoCrMo Alloy (ISO 5832-4)/ CoCrMo Alloy (ISO 5832-4)
 - Implant Steel (ISO 5832-1)/Titanium Alloy (ISO 5832-3)
 - Implant Steel (ISO 5832-1)/Implant Steel (ISO 5832-1)
- Components that are incompatible with each other:
 - Implant Steel (ISO 5832-1)/ CoCrMo Alloy (ISO 5832-4)
5. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
6. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning, and must inform the patient.
7. These products are for single use only. Do not re-use the product.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
2. The surgeon must warn the patient that the prosthesis cannot replace the normal healthy bone, prosthesis can be broken or damaged as a result of some movements or trauma, the expected product life is limited, and replacement can be needed in future.
3. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.
4. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.
5. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.
6. Incorrect selection or positioning of components can result in insufficient femoral range of motion because of femoral squeeze and peri-articular calcification.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

7. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

8. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred. For this purpose, fluoroscopy can be required.
2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.
3. Correct section of the shell and stem based on the length of the neck is very important. Muscular laxity and/or incorrect placement of components can lead to loosening, partial displacement and/or break of implants. Longer neck length and placement with inward inclination will increase the tension that should be borne by the stem. Components must be tightly implanted in place using surgical hand tools.
4. As it is not possible to see the internal tension of a prosthesis implanted before, bending or fractures in the early period can be seen; therefore, prosthesis used before must never be re-used.
5. During the implantation, femoral or acetabular perforation/fractures can occur because of trauma or overloading, particularly in the presence of poor bone structure.
6. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.
2. The patient must be warned that the prosthesis cannot replace the normal health bone, the prosthesis can be broken or damaged as a result of some movements or trauma, the expected product life is limited, and replacement can be required in future.
3. Abrasion of the components can result in the increase of loosening and damaging of the bone.
4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending or cracking of components.
5. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.
6. Recommendations of surgeons in the postoperative period are of utmost importance.
7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarct can occur.

ZİMED MEDİKAL
INSTRUCTIONS FOR USE OF REVISION HIP PROSTHESIS
Submitted to the Attention of Surgeons
(EN)

Document No: IFU.10/00
14.06.2022

2. Damaging of blood vessels or hematomas can be seen.
3. Delayed wound healing: Deep wound infections (early or late) that might require the removal of the prosthesis can occur. Arthrodesis of the relevant joint or amputation of the extremity can be required in rare cases.
4. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
5. Unwanted lengthening or shortening of the extremity can occur.
6. Fatigue fractures of prosthesis components can arise from trauma, exhausting physical activity, misalignment, defective implantation, service period, loss of fixation, nonunion or overloading.
7. Trochanteric nonunion can result from insufficient re-connection and/or load-bearing in the early period. Prosthesis can collapse.
8. Trochanteric avulsion can result from excessive muscular tension, weight loss in the early period or inadvertent intra-operative weakening.
9. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
10. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
11. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
12. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body can change according to the daily activity of the patient. Product life of this implant has been determined as 15 years. Implants can be removed according to the surgeon's decision.

STERILIZATION:

Implants are submitted as sterile. Sterilization method will be shown on the packaging. Minimum 25 kGy gamma irradiation is applied to all the implants sterilized with gamma irradiation. Sterile products will remain so unless the integrity of the package is broken. The products are presented as sterile. **Re-sterilization** of the product is not recommended. Please inform the manufacturer about any products with impaired sterility although offered as sterile.

SHELF LIFE:

The shelf life for **Revision Hip Prostheses** sterilized with Gamma sterilization is 5 years. Please check the expiry date. Please contact the manufacturer for products expired. Non-sterile products do not have shelf lives.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: *Description of Symbol*)

SEMBOL AÇIKLAMALARI	
	Logo of the Firm
	"1984" Notified Body Number
	Medical Device
	Please Read the Instructions for Use
	Lot Number
	Product Reference/Catalog Number
	UDI Number (Single Device Identification Number)
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation.
	Sterilized with ethylene oxide
	Used only with prescription.
	Implants are for single use.
	Do not re-sterilize.
	Do not use damaged packages
	Keep away from direct sunlight
	Keep away from precipitation.

	Fragile: Handle with care
	Temperature limits
	Caution

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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