www.zimed.com.tr

zimed®

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RU

ZİMED MEDİKAL

INSTRUCTIONS FOR USE OF SPINAL SYSTEMS

Submitted to the attention of Surgeons (EN)

Document No: IFU.03/00

14.06.2022

(EN) PRODUCT DEFINITION:

Spinal Systems are designed to be used in acute and chronic instabilities, deformations of thoracic, lumbar and sacral vertebral areas for fusion treatment, in patients with mature skeletal systems, for the fixation of skeletal structure and to ensure the stabilization of spinal segments, after failed fusion treatments, and in fractures, dislocations, scoliosis, kyphosis and in deformations of spondylolisthesis. Spinal system products are dispatched as non-sterile.

Polyaxial Screw is made from Ti 6Al4V ELI ASTM F136 (5832-3) material with the purpose of using in reduction of vertebral fractures, displacements or traumas and to ensure cervical fusion stabilization, and has a unique thread structure. Polyaxial screws can be positioned with 60° angular movement.

Monoaxial Screw is made from Ti 6Al4V ELI ASTM F136 (5832-3) material with the purpose of using in reduction of vertebral fractures, displacements or traumas and to ensure cervical fusion stabilization, and has a unique thread structure. Monoaxial screw do not have the capability of angular movement.

Spondylolisthesis Polyaxial Screw is the screw capable of angular movement, and is made from Ti 6Al4V ELI ASTM F136 (5832-3) material to align the displacements of vertebral discs.

Spinal Rod is made from Ti 6Al4V ELI ASTM F136 (5832-3) material to ensure stabilization of screws and to support the spine in the treatment of scoliosis, trauma and in fusion treatment.

Transverse Connector Crosslink is made from Ti 6Al4V ELI ASTM F136 (5832-3) material, and ensures stabilization by creating the connection between rods.

Transverse Connector Hook is made from Ti 6Al4V ELI ASTM F136 (5832-3) material in one single type with standard measures, and ensures the connection between the spinal rod and transverse connector crosslink.

Angular Transverse Connector is made from Ti 6Al4V ELI ASTM F136 (5832-3) material to ensure the connection between the spinal rod and the transverse connector crosslink, and is capable of angular positioning.

Please refer to the product catalog for the reference numbers and names of products.

INTENDED USE AND INDICATIONS:

Spinal Systems are used in patients with mature skeletal structures in order to assist bone healing by ensuring stabilization and preventing mobility of spinal segments, and in the treatment of acute and chronic thoracic, lumbar and sacral vertebral injuries and instabilities.

In addition, this device is used in severe cases of spondylolisthesis of L5-S1 vertebrae (grade 3 or 4) in patients with mature skeletal systems, and they are removed following the completion of spontaneously-occurring bone healing of implants in sacral spine. In such indications, level of screw fixation in pedicles include only the vertebrae between L3 and S1.

Spinal Systems are for pedicle and non-pedicle fixation, and is designed only for use in the posterior side. When used in patients with mature skeletal systems for the fixation of pedicles in the non-cervical posterior spine, the system is indicated for one or more of the following:

Degeneration of thoracic and lumbar discs,

Spondylolisthesis

Trauma (such as fractures or luxation/dislocation)

Deformities (scoliosis, kyphosis and/or lordosis, Scheuermann disease)

Tumors

Spinal stenosis

Failed fusion revision

Pseudoarthrosis

Iatrogenic or primary instabilities

In addition, this system indicated for the treatment of advanced spondylolisthesis between the autogenous bone graft and the L5-S1 joint (grade 3 or 4) and treatment of fusion treatment in lumbar and sacral vertebrae and removal from the healed tissue.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

CONTRA-INDICATIONS:

Active infections, manifest or occult infection or manifest local inflammation in or around the affected site, existing fracture fixation device implanted before, hypovolemia, hypothermia, coagulopathy, obesity, mental disorders preventing cooperation with the rehabilitation

regime, advanced osteoporosis, poor bone quality or amount, disorders in the medullar canal, poor blood circulation, and previous infections will delay healing.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products of different manufacturers must not be used together.
4. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
5. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning, and must inform the patient.
6. These products are for single use only. Do not re-use the product.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
2. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.
3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.
4. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.
5. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.
6. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative:

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.
2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.
3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.
4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and

protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. The patient must be warned that the implant cannot replace a normal and healthy bone, the implant can break or be damaged as a result of some movements and traumas, its expected life is limited, and replacement can be necessary in the future.
3. Abrasion of the components can result in the increase of loosening and damaging of the bone.
4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending or cracking of components.
5. Fatigue fractures of implants can arise from trauma, exhausting physical activities, wrong alignment, defective placement of the implant, service time, loss of fixation, nonunion or overloading.
6. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.
7. Recommendations of surgeons in the postoperative period are of utmost importance.
8. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Damaging of blood vessels or hematomas can be seen.
2. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
3. Unwanted lengthening or shortening of the extremity can occur.
4. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
5. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
6. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
7. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body is 12 to 18 months, and must be removed according to the doctor's decision.

STERILIZATION:

Implants Put on the Market as Sterile: Spinal System products are submitted as sterile or non-sterile. Sterile products will remain sterile until the package integrity is broken. There are two types of sterilization as sterilization with gamma irradiation or ethylene oxide. Sterilization method will be shown on the packaging. Minimum 25 kGy gamma irradiation is applied to all the implants sterilized with gamma irradiation.

Please inform the manufacturer about any products with impaired sterility although offered as sterile.

Implants Put on the Market as Non-Sterile: Products put on the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in.

All the Spinal Systems products and surgical hand tools used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

Autoclave Sterilization: Non-sterile products put on the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel : 0 342 238 43 44 Fax : 0 342 238 44 11 Web:www.zimed.com.tr E-mail: info@zimed.com.tr

ZİMED MEDİKAL
INSTRUCTIONS FOR USE OF SPINAL SYSTEMS
Submitted to the attention of Surgeons
(EN)

Document No: IFU.03/00
14.06.2022

AUTOCLAVE STERILIZATION	
Parameters	Autoclave with Pre-vacuum
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying Time	15 minutes

Note: Exposure time is the period to be applied starting from the moment that the autoclave reaches the exposure temperature.

The following table shows the autoclave types as compared according to certain sterilization parameters.

PRODUCT TYPE	RECOMMENDED STERILIZATION METHOD
Spinal Systems	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vakum Otoklav): <i>If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)</i>

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

SHELF LIFE:

The shelf life for Spinal Systems sterilized with Gamma sterilization is 5 years, while the same for those sterilized with ethylene oxide is 3 years.

Please check the expiry date. Please contact the producer for products expired.
Non-sterile products do not have shelf lives.

DISPATCH OF THE PRODUCTS:

The Spinal Systems presented as sterile or non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the producer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: Description of Symbol)

SYMBOL DESCRIPTIONS	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single-Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation.
	Sterilized with ethylene oxide
	Used only with prescription.
	Implants are for single use
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight

	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution
	Non sterile

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel : 0 342 238 43 44 Fax : 0 342 238 44 11 Web:www.zimed.com.tr E-mail: info@zimed.com.tr

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