



Trauma | Arthroplasty | Sports Medicine | Bone Substitutes

## PFNA SURGICAL TECHNIQUE

## CONTENTS

### **3 Lag Screw Specification**

#### 4 Fracture

##### Entry Point

#### 5 Approach

### **6 Opening the Proximal Femur**

##### Kirschner Wire

##### Proximal Drill Guide

#### 7 Proximal Reamer

### **8 AWL For First Entry**

##### Optional Option

##### Guide Wire

#### 9 Optional Option

##### Proximal Reamer

### **10 Nail Insert**

##### Preparing the nail

##### Guide wire

### **11 Inserting Nail**

##### Control

### **12 Connection Devices**

##### Carbon Kishner Guide Optional Device

### **13 Preparing Lag screw**

##### Lag Screw Guide Point Marker

### **14 Lag Screw Guide - Kirschner Guide**

### **15 Prepare drill**

##### Drilling

### **16 Lag Screw**

##### Prepare Lag Screw

##### Lag Screw insert

### **17 Compression**

##### Lag Compression Screw

### **18 Preparing lag end cup**

##### Attaching lag end cup

### **19 Distal Screw**

##### Removing targeting devices

##### Distal Targeting Device Attachment

### **20 Removing targeting devices**

### **21 Nail End Cap**

##### Removing handle

##### Nail End cap attachment

### **22 Interclaw Lag Screw with Blade Proximal Femoral Nail**

### **23 Interclaw Lag Screw**

##### Preparing The Lag

##### Lag Screw attachment

### **24 Blade oppening preparing**

### **25 Compression**

##### Lag compression screw

### **26 Lag end cap**

### **28 Finish**

### **29 Extract**

##### Nail extracting devices

##### End Cap Removing

### **30 Nail extracting preparing**

##### Lag screw removing (Lag end cap)

##### Lag screw removing (Preparing for blade close)

### **31 Blade closing**

### **32 Distal screw removing**

### **33 Anti-Rotation Proximal Femoral Nail**

#### 34 Guides

##### Kirshner Wire

#### 35 Measurement

##### Anti-Rotation device

### **36 Drilling&Simple Lock Lag Screw**

### **37 Compression Screw**

### **38 Compression**

### **39 Appendix**

#### 40 Long Nail, Determining the length of the nail & thickness

##### Devices for first entry

#### 41 First entry & drilling

#### 42 Reamer

#### 43 Determine Nail Lenght

#### 44 Disinfection

## Lag Screw Proximal Femoral Nail



### PFN-A Intramedullary Nail

The PFN-A Nail is used in the treatment of unstable intertrochanteric fractures. Intertrochanteric femur fractures are common in the elderly population. This is because of osteoporosis. Due to the decrease in bone quality and deterioration of its microstructure, fractures often develop with very low-energy trauma

Lag is the most important feature of the nail. There are three types.

They come to the forefront with their features such as screw thread structures, compression, blade features and they have their own unique application forms.

### Lag Screw

Can be used in low-energy unstable intertrochanteric fractures and in patients with non osteoporosis and younger

### PFN-A Nail Technical Specifications

- Proximal diameter Ø16mm
- Distal diameter Ø10mm - Ø11mm - Ø12mm - Ø13mm - Ø14mm
- Proximal – Distal angle 5°
- Lag screw center angle 55° to proximal body
- Cannula diameter Ø3,7mm for all diameters of PFN-A nails
- Distal antirotation (Ø2,5mm X 25mm)
- Ø5mm locking screw dynamic locking screw

REF. NO	LENGTH
4582-0080	80
4582-0085	85
4582-0090	90
4582-0095	95
4582-0100	100
4582-0105	105
4582-0110	110
4582-0115	115
4582-0120	120


## Fracture

Evan's classification system divides intertrochanteric fractures into stable and unstable fracture patterns. The distinction between stable and unstable fractures is based on the integrity of the posterior medial cortex. Other intertrochanteric fracture classifications are variations of Evan's classification, including AO.

In general, when the posterior medial cortex is fragmented, fractures are considered unstable due to the possibility of the fracture collapsing into varus and retroversion.(Fig.1)


It is preferred in simple and complex intertrochanteric fractures.

Appropriate lag screw can be selected according to fracture severity and bone quality.



## Entry Point


The design of the nail is designed according to the slope. It will give the appropriate angle when placed from the right point.






## Approach

A 3-5 cm skin incision is made a few centimeters proximal to the tip of the greater trochanter. It is located in the proximal extension of the curved axis of the femoral shaft.



Pay attention to the systems that need attention



## Opening the Proximal Femur Kirschner Wire

Send the wire from the entry point you have determined with the motor.



## Proximal Drill Guide

Place the guide over the kirschner wire



Proximal Drill Guide Handle  
Below the carving process to preserve  
the tissues.  
Ø17mm  
Ref:9455-0051



## Opening the Proximal Femur Proximal Reamer

Insert the proximal reamer  
into the guide  
Remaining up to the stopper on  
the reamer



Proximal Reamer  
Ø17mm Cannulated Reamer  
Ref:9455-0047



Remove proximal reamer, K  
wire, Guide





## AWL For First Entry Optional Option

Open the cortex with awl



AWL  
Ø3.5mm  
Ref:9455-0016



## Guide Wire

Find the canal with  
guide wire



Gripper



## AWL For First Entry Optional Option

Remove Awl  
Guide wire is stands in canal



---

## Proximal Reamer

Insert the proximal reamer  
into the guide over guide wire  
Remaining up to the stopper on  
the reamer



## Nail Insert

### Preparing the nail

For the nail, the fracture is reduced and the channel is opened. The length of the nail was also determined. To be able to insert the nail into the canal, it must be combined with the holder. The connecting screw is attached to the holder. It is fixed with screwdriver



Targeting Device Handle  
Ref:9455-0030



T Holder 8 Allen  
Ref:9455-0018




Cannulated Targeting  
Connection Screw  
M10xØ4.6x42.5mm  
Ref:9455-0041




## Guide wire

Attach it to the gripper with guide wire. Insert the guide wire through the opened inlet channel, and advance it through the opened channel.



## Inserting Nail

After it is attached to the holder the driving apparatus and the hammer swing shaft are combined. It's advanced in the channel with the help of hammer



## Control

Check by inserting the Kirschner wire into the holes on the guide. At this stage, you will also determine the size of the end cap that you will insert in the following stages.






## Connection Devices

Attached insertion Handle  
with Proximal Lag Screw  
Targeting device handle





## Carbon Kishner Guide Optional Device

This optional product allows you to check in advance the length of the Lag screw and whether it will be in the right place with the Kirshner wire you will attach.



## Carbon Kishner Guide Optional Device



## Preparing Lag screw Lag Screw Guide Point Marker

Insert the group sleeve and mark with the point maker and remove it

Lag Screw Drill Sleeve  
Ref:9455-0024

Point Marker  
Ref:9455-0026



Kirschner Sleeve  
Ref:9455-0025



## Preparing Lag screw Lag Screw Guide Kirschner Guide

Insert the Kirshner wire with motor

Remove Kirshner Sleeve  
Control under the image  
place of Kirchner wire. Right  
place important Beacuse Lag  
screw will send over this wire.




## Lag Screw Guide Kirschner Guide

Use the Lag Screw gauge  
determine screw lenght over  
the Krishner wire




Lag Screw Depth Guide  
Ref:9455-0046




## Preparing Lag screw

### Prepare drill

Prepare Ø10.8xØ3.2mm Cannulated Drill Bit) for drilling. Adjust ( drill stoper) according to prestage size measurement



#### NOTE



## Drilling.

Insert the Cannulated drill bit into the Drill sleeve and make drill to the stopper






## Lag Screw Prepare Lag Screw



Lag Screwdriver  
Ref:9455-0000



### Lag Screw Lengths

80  
85  
90  
95  
100  
105  
110  
115  
120



Selecting the Right Screw Length  
Based on the Measurement during the  
Drilling Phase



"The screw that fits the  
determined size is attached  
to the screwdriver."

Lag Screw Driver  
Ref:9455-0000

## Lag Screw Lag Screw insert



The guide is used to insert the  
screw with screwdriver

Red Marker is used for 0 point  
its show of where the end of  
screwdriver



## Lag Screw Compression


Compression is applied as needed based on the requirements.



## Lag Screw Lag Compression Screw



Combine the compression screwdriver with the features of the compression screw. In this way, the screw is prevented from falling out during the installation phase.



Send the lag lag locking screw through handle to proximal of the nail




## Lag Screw Preparing lag end cup


Before Remove the lag screwdriver and Kirschner wire.



use allen key and attach lag end cap



## Attaching lag end cup



send the lag end cap through to the guide





## Distal Screw Removing targeting devices

The lag screw installation process is completed. Sleeve and aiming device are removable



## Distal Targeting Device Attachment



Dynamic-Static Locking  
Distal Target Arm  
Ref:9455-0055




Ø4.0mm Sleeve  
Ref:9455-0019

Ø4.0mm Guide  
Ref:9455-0020

ØØ4.0mm Point Marke  
Ref:9455-0021



Ø4.3 x 300 Drill  
Ref:9455-0008



Depth Guide  
Ref:9455-0055

For the distal screw, a distal screwing guide is attached. It is fixed with a connecting screw. Secure with L allen



Static: its used for 170to200mm lenght short nail

Static: its used for 240mm lenght short nail

Dynamic: its used for 170to200mm lenght short nail

Dynamic: its used for 240mm lenght short nail

Decide for dynamic or distal locking

## Distal Screw Removing targeting devices

If necessary, a dynamically static distal screw can be sent. Marking is done and marker is removed. Drilling with motor and remove sleeve decide screw length with depth guide. Send screw



Send Sleeve



Point marked and remove



Drilling




Remove Sleeve



Decide screw length with depth guide



Send the Screw



The distal screw head is grooved and can be held with a screwdriver.



Remove sleeve and distal targeting device

## Nail End Cap Removing handle




T Holder 8 Allen  
Ref:9455-0018

Remove (distal targeting device & handle device) before sending end cap



## Nail End cap attachment




REF. NO	LENGTH
4562-0000	0
4562-0005	5
4562-0010	10




Anti-Rotation Screwdriver  
Ref:9455-2065

T Allen Key  
Ref:9456-0003



There are 3 types of end cups. Install the appropriate one.



Complete the process by inserting the End Cap screw with a screwdriver. Make the final check under Imaging and complete the process.

## Interclaw Lag Screw with Blade Proximal Femoral Nail



### PFN-A Intramedullary

The PFN-A Nail is used in the treatment of unstable intertrochanteric fractures. Intertrochanteric femur fractures are common in the elderly population. This is because of osteoporosis. Due to the decrease in bone quality and deterioration of its microstructure, fractures often develop with very low-energy trauma

### Interclaw Lag Screw with Blade

Can be used in high-energy unstable intertrochanteric fractures and in patients with osteoporosis and older

### PFN-A Nail Technical Specifications

- Proximal diameter Ø16mm
- Distal diameter Ø10mm - Ø11mm - Ø12mm - Ø13mm - Ø14mm
- Proximal – Distal angle 5°
- Lag screw center angle 55° to proximal body
- Cannula diameter Ø3,7mm for all diameters of PFN-A nails
- Distal antirotation (Ø2,5mm X 25mm)
- Ø5mm locking screw dynamic locking screw



## Interclaw Lag Screw

REF. NO	LENGTH
4592-0080	80
4592-0085	85
4592-0090	90
4592-0095	95
4592-0100	100
4592-0105	105
4592-0110	110
4592-0115	115
4592-0120	120



## Interclaw Lag Screw Preparing The Lag.

The blade, which comes mounted inside the lag screw, ensures its retention when opened.



Note: The stages up to the bladed lag screw are applied with the same procedures as the lag screw.

Prepare lag screw driver for blade and lag screw



Lag Screw Driver  
Ref:9455-0000


## Lag Screw attachment

Send Interclaw lag screw with screw driver over the Kirschner Wire and remove wire



## Interclaw Lag Screw Blade opening preparing.


The blade, which comes mounted inside the lag screw, ensures its retention when opened.



## Blade opening preparing.



Open Blade



## Interclaw Lag Screw Compression

If you need compression, perform compression by turning it with an allen wrench, and check under visualization.



## Lag compression screw



Screwdriver  
4.0mm  
Ref:9455-0055

Prepare lag locking screw with screw driver and attach

it is used for stabilization of the lag to prevent it from moving in any situation





## Interclaw Lag Screw Lag end cap



Prepare lag end cap and send it with screwdriver



Lag end cap  
Ref:9455-0055



Remove the screwdriver



## Nail End Cap Lag end cap


To install the nail end cap, remove the c arm .



T Holder 8 Allen  
Ref:9455-0018



REF. NO	LENGTH
4562-0000	0
4562-0005	5
4562-0010	10




Prepare nail end cap

and then install the nail end cap


**Note:**  
Distal screw attachment procedure same for all short nail technique. Look for distal screw stage

finish and provide necessary controls



## Finish

Complete the procedures by providing the necessary checks and complete the surgery.



The procedure is complete. Depending on the patient's condition and the doctor's decision, the implants can be removed when healing occurs. For this procedure, you can review the removal procedures section.

## Extract Nail extracting devices


After the necessary checks, the nail may need to be removed due to the patient's request and the doctor's approval, and similar reasons. This section shows the removal of the nail and the removal of the bladed lag screw.

The entire set is not needed to remove the nail. Some of the tools in the set are used. Depending on the situation, only the tools used for removal are sent to the surgery.


Device for  
extracting



## End Cap Removing




Remove Nail end cap with Anti-rotation screw driver



Remove lag locking screw with screwdriver 4.0mm


## Extract Nail extracting preparing.

In order to remove the nail, insert the rod at this stage (before removing the lag screw). In this way, you will prevent the nail from falling into the channel when you remove the lag screw.



## lag screw removing (Lag end cap)

Lag end cap must be removed before removing lag screw




## lag screw removing (Preparing for blade close)



## Extract Blade closing.

Attach lag screw driver to lag screw and others



Rotate handle and blade will close.

You can remove lag screw securely now.





## Extract Distal screw removing.

Remove distal screw with screw driver



Remove nail with moving hammer over the rod. Nail is now removed. Complete other procedures and check-ups required to complete the surgery.





## Anti-Rotation Proximal Femoral Nail



### PFN-A Intramedullary

The PFN-A Nail is used in the treatment of unstable intertrochanteric fractures. Intertrochanteric femur fractures are common in the elderly population. This is because of osteoporosis. Due to the decrease in bone quality and deterioration of its microstructure, fractures often develop with very low-energy trauma

### Simplelock Lag Screw

may be preferred in high-energy and unstable fractures if more compression is needed

### PFN-A Nail Technical Specifications

- Proximal diameter Ø16mm
- Distal diameter Ø10mm - Ø11mm - Ø12mm - Ø13mm - Ø14mm
- Proximal – Distal angle 5°
- Lag screw center angle 55° to proximal body
- Cannula diameter Ø3,7mm for all diameters of PFN-A nails
- Distal antirotation (Ø2,5mm X 25mm)
- Ø5mm locking screw dynamic locking screw

REF. NO	LENGTH
4792-0060	60
4792-0065	65
4792-0070	70
4792-0075	75
4792-0080	80
4792-0085	85
4792-0090	90
4792-0095	95
4792-0100	100
4792-0105	105
4792-0110	110
4792-0115	115


## Anti-Rotation Proximal Femoral Nail Guides

Decide Lag screw right place with guide

Attach lag drill guide through inside  
targeting device and Kirshner guide

### Kirshner Wire

Send the Kirshner Wire with motor



Kirshner Wire Two Type  
Curved & Normal  
Threaded

## Anti-Rotation Proximal Femoral Nail

### Measurement

Send the Measurement Guide over the  
Kirshner Wire



Send the Measurement Guide over the  
Kirshner Wire



Drill with Lag drill for compression screw,  
Drill -5



### Anti-Rotation device

Attach anti rotation device




Anti rotation device  
ref:



## Anti-Rotation Proximal Femoral Nail

### Drilling & Simple Lock Lag Screw

Remove Kirshner guide




Adjust the stopper according to the measurement and perform drilling.

Select and install the lag screw according to the measurement length



Choose  
110 mm



## Anti-Rotation Proximal Femoral Nail Compression Screw

Send the screw until the number zero.



Remove antirotation device



Choose a 105 mm screw. Because during the drilling phase, we drilled 5 mm less than the 110 mm we measured.



**Choose  
105 mm**






## Anti-Rotation Proximal Femoral Nail

### Compression

Make sure that the yellow mark on the compression screw is inserted until it matches the zero number on the lag screw.


If you are not going to tighten, you can remove the screwdrivers and move on to the next steps. If you are going to tighten, turn the lag screwdriver to 5. If you need more tightening, you can turn it up to 10.

Number 0  
no need compression




compression example


Number 0  
need compression



Number 5  
Try number 5 but  
insufficient compression



sufficient compression was  
achieved



## Appendix

Long Nail, Determining the length of the nail & thickness




The long Pfna nail is used in the treatment of intertrochanteric fractures as well as fractures occurring in the shaft.

## Appendix


Long Nail, Determining the length of the nail & thickness

### Devices for first entry



Gripper  
Ref:9413-0002

Guide Wire  
02.5/03.3x1000mm



AWL  
Ref:9455-0016



## Appendix

Long Nail, Determining the length of the nail & thickness

### First entry & Drilling

Open the cortex with awl

Find the canal with  
guide wire





Remove Awl



use proximal reamer sleeve and make drill  
with proximal reamer

remove sleeve






## Appendix

Long Nail, Determining the length of the nail & thickness

### Reamer



**Flexible reamer**  
Ref:9426-0008



**Reamer Heads**  
Ref:9455-0033-0039



Combine reamer and start with small size





## Appendix

Long Nail, Determining the length of the nail & thickness

### Determine Nail Length

You can determine the thickness of the screw with reamer tips. Start from the smallest size and carve towards the largest size. In this process, you can decide the thickness of the nail



Short PFNA

Long PFNA

The appropriate nail size can be determined from the lines on the reamer.

### Short

10x170	12x200
10x200	12x240
10x240	13x170
11x170	13x200
11x200	13x240
11x240	14x170
12x170	14x200
	14x240



### Long

Left	Right
10x340	10x340
10x360	10x360
10x380	10x380
11x340	11x340
11x360	11x360
11x380	11x380
12x340	12x340
12x360	12x360
12x380	12x380
13x340	13x340
13x360	13x360
13x380	13x380

## Disinfection

### DEVICE CLEANING CONDITIONS

Do not use metal brushes or rubbing pads during Decontamination of the tools should be performed immediately after the surgical procedure is completed. Contaminated tools must not be allowed to dry before reprocessing.

Excessive blood or debris must be removed in order to prevent the drying on the surface. All users must be qualified staff with documented evidence of training and competence. Training should include the current guidelines, standards and hospital policies. Even if they are made of high-grade stainless steel, the surgical tools must be thoroughly dried in order to prevent rust formation. Prior to sterilization, all the tools should be examined for the cleanliness of the lumens of the joints of the surfaces. manual cleaning process. Use cleaning agents with low-foam surfactant to be able to see the tools in the cleaning solution. Rinse the cleaning materials easily from the tool in order to prevent residue formation.

Mineral oil or silicon lubricants should not be used on Zimed tools. Neutral pH enzymatic and cleaning materials are recommended for cleaning the reusable instruments. It is very important to neutralize and rinse the alkaline cleaning materials thoroughly from the tools. Anodized aluminum should not contact with certain cleaning or disinfectant solutions. Avoid strong alkaline cleaners and disinfectants and solutions containing iodine, chlorine or certain metal salts.

### Manual Cleaning/Disinfection

Prepare the enzymatic and cleaning materials at the dilution rates and temperatures as recommended by the manufacturer. New solutions should be prepared when the existing solutions are heavily contaminated. Place the tools in the enzymatic solution so that they are completely immersed. Operate all the movable parts so that the detergent contacts with all the surfaces.

Keep in the fluid for minimum 20 min. Use a nylon, soft-bristled brush to gently rub the tools until all visible debris is cleaned. Pay particular attention to the accessible areas and use a suitable bottle brush. In order to remove the dirt in the open springs, coils or flexible parts, wash the recesses with plenty of cleaning solution. Rub the surface with a scrubbing brush to remove all the visible dirt from the surface and the recesses. To ensure that all the recesses are cleaned, turn the component while rubbing. Remove the tools and rinse them for minimum 3 min. under running water. Pay particular attention to the cannulas and use a syringe to pass the fluid through the hard-to-reach areas. Place all the tools that are completely immersed in water, in an ultrasonic unit containing the cleaning solution. Operate all the movable parts so that the detergent contacts with all the surfaces. Expose the tools to sonification process for minimum 10 min..

Remove the tools and rinse with deionized water for at least 3 minutes or unless all the blood or dirt traces are eliminated in the rinsing water. Examine the tools under normal light to verify that visible dirt is removed. If visible dirt is present, repeat the above mentioned sonification procedure and the rinsing steps. Remove the excessive moisture on the tool with a clean, absorbent, lint-free cloth.

### Combination Manual / Automated Cleaning and Disinfection

Prepare the enzymatic and cleaning materials at the dilution rates and temperatures as recommended by the manufacturer. New solutions should be prepared when the existing solutions are heavily contaminated. Place the tools in the enzymatic solution so that they are completely immersed. Operate all the movable parts so that the detergent contacts with all the surfaces. Keep in the fluid for minimum 10 min. Use a nylon, soft-bristled brush to gently rub the tools until all visible debris is cleaned. Pay particular attention to the accessible areas and use a suitable bottle brush. A sonicator will help to clean the instruments thoroughly. The use of a syringe or a water fountain will facilitate passing of the liquid from the low-spaced areas and difficult-to-access areas. Remove the tools from the enzyme solution and rinse them for minimum 1 min. under deionized water. Place the tools in a suitable washer / disinfectant basket and perform a standard washer / disinfectant cycle. Specific minimum parameters are essential for a complete cleaning and disinfection. These parameters are given in a below mentioned table.

### Automated Cleaning and Disinfection

Automated washing / drying systems are not recommended as the only cleaning method for surgical tools. An automated system can be used as a follow-up operation after manual cleaning. To ensure an effective cleaning, tools must be thoroughly examined before sterilization. For detailed information on Washing and Disinfection see

### Specific minimum parameters used for a complete cleaning and disinfection:

Definition	
1	Pre-washing for 2 minutes with cold tap water
2	enzyme spray for 20 seconds with hot tap water
3	Immersion in enzyme after 1 minute
4	rinsing for 15 seconds with cold tap water (Should be repeated twice)
5	Washing with detergent for 2 minutes with hot tap water
6	rinsing for 15 seconds with hot tap water
7	Rinsing with 10 seconds with optional lubricated purified water
8	Drying for 7 minutes with hot air

Note: Follow the instruction of the washer/disinfectant manufacturer

*Zimed Medical, as the manufacturer of this device, and their surgical consultants do not recommend this or any other surgical technique for use on a specific patient. The surgeon who performs any implant procedure is responsible for determining and utilizing the appropriate techniques for implanting the device in each individual patient. Zimed and their surgical consultants are not responsible for selection of the appropriate surgical technique to be utilized for an individual patient.*



ISO 9001  
ISO 13485 CE 1984